

Product datasheet for **MC211316**

Tmem158 (NM_001002267) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tmem158 (NM_001002267) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Tmem158
Synonyms:	2310037P21Rik; Ris1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC211316 representing NM_001002267 Red =Cloning site Blue =ORF

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGCCCTTCTCGCCGCGCTACTGGCCGCCCTGCCAGCTGCCGCCGGCGCACGGCGGGGCCACGG
ATGCGCCCGCCTTGCAGGGACGCCCCAAACGCCTCGGCCAATGCGTCTTACCAACGAACACTCCAC
CCCGCTGCTGGCATCGGCGGCATCGGCGCCCCAGAGCGCGGGGCCAGAGGAGGCTCCGGCGGCG
CCCTGCAACATCAGCGTCAACGACAGATGCTGAGCTCGTCTGCTGCGCTGGGGCCGCCCGGGGCT
TGCAGTGGACCTGCTGCTTCTCTACCAACGCGCACGGCCGCGCTTCTTCGCCCGCCCTTCCACCG
CGTCGGGCCACCACTGCTCATCGAGCATCTGGGGCTGGCGGGGGCGGCGCAGCAAGACTGCGCCTG
TGGCTGGGCTGCGGCTGGGTGCGCGGCCCTGCGGGCTCCCGCCGGGGTCCCACCGCGCTGCCTGCCT
ACCCGGCGGCCGAGCCGGGGCCGCTGTGGCTGCAGGGCGAGCCTCGTCACTTCTGCTGCCTGGACTTCAG
CCTGGAGGAGCTGCAGGGCGAGCCGGGCTGGCGGCTGAACCGCAAGCCATCGAGTCTACGCTGGTGGCC
TGCTTCATGACCCTGGTTCATCGTGGTGTGGAGCGTGGCCGCCCTCATCTGGCCGGTGCCCATATTGCCG
GCTTCTGCCAACGGCATGGAGCAGCGCCGACCACCGCCGGCGCCCCGCGGCCGCTCCTGCAGCCGT
ACCCGCGGGGACTACCGCCGCTGCCGCCCGCTGCTGCAGCTGCTGCTGCTGCTGCCCGGCTGTCAAC
TCGGGGTGGCTCCTAAGTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGTTTAA

Chromatograms:	https://cdn.origene.com/chromatograms/ja1654_a05.zip
Restriction Sites:	Sgfl-MluI
ACCN:	NM_001002267



[View online »](#)

Insert Size:	861 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC139427</u> , <u>AAI39428</u>
RefSeq Size:	969 bp
RefSeq ORF:	861 bp
Locus ID:	72309
Cytogenetics:	9 F4