

Product datasheet for MC211347

Pih1d2 (NM_028300) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids

Tag: Tag Free

Symbol: Pih1d2

Synonyms: 2700059L22Rik

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC211347 representing NM_028300
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGACGGTGTCTCAAAGGGTCTTCTTACCCACATTTCCAGTTCTGGAATATGCTCGATGATTTGGCTG
AAAATGACCCCGAGCGCTACCGGAACCTTCATTCAGCAGGAGTTGAAAGATGGCAAACAGCTCTGTGTCAA
CCCAGAACCTCAACTCTGCATACAGACCAAGATCCTAAAACCAAACGAAAAAGTACTTTTATCAACCTA
TGTCATGGGAAAGGATTCCCGCTCCGAGTCAGCCACTCGTCTGTGCCTGTAAGTGTGGCAGACCAG
AAGACTCCGCTGAGGCATCAGATGCTTATACCATCATTGATGTTGCCTACAATCCTGGTGTCTGCAAGC
GGCAGAAAAAGACCAAGGGATCAAAGATCAGTTAATAAGAATGGCCATGCTGTGCATCGAGGAGCGACTC
CAGTTCACACTTGACATTCTACCATTTAACTAGTTTACTGAAAGGAAGCATCCAGAGAATGAAAG
AAAGTCTGATGGGCATTAAACTGACTTCCAGACCTAAAAGCGATGACGAGAACAGAAAACTCTTGC
GAGGATCAGAAGCAGCACTGTGAGTGAACCTAATCACCTTCTGAGGTGCTATTGACAAAAAGCAAGCG
TCGGCCAAAGGACGATGTCTTATAGAAGAGATTCTAGTTCGGAATCCAGGTGGAGTTAAGAAACCTG
CGTACGAATTAAGGTTGTAAAGGATCGGAACGAGAAACCTCTGAAATGAGCTGAAAGTCGAGTTACC
CGGGATTAAGTCAGTCTCTCTGTGAACCTAGTGTCTGAGGTTGATATATTGATTGAGGTCTCTGAG
ATGTACCGATTATGCCTGAATCTTCAGAATCGATTAATACAGAAATGACTACAGCAAAATTTGTCAAAA
ATAAATCTGCATTAATCATTACAATGCCCTGGCATAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_028300



Insert Size:	948 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_028300.2</u> , <u>NP_082576.1</u>
RefSeq Size:	1268 bp
RefSeq ORF:	948 bp
Locus ID:	72614
UniProt ID:	<u>Q8CHR9</u>
Cytogenetics:	9 A5.3