

## Product datasheet for MC211391

### 3110079O15Rik (NM\_028473) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	3110079O15Rik (NM_028473) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	3110079O15Rik
Synonyms:	Snorc
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC211391 representing NM_028473 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCATCTTGTCTGGCCCTGCGCTGGCGTGCTGCTCATCTCTGGGTCTGGCCCTGCGGTGCTCA  
CAGCAGAGGGTCCCCAGGAGCCGATCCCACCCTGTGGAACGAGCCTATTGAATTGCCCTCTGGTGAAGG  
TCCCTGGAGAGCACAAGCCACAACCAGGAATTTGCTGTCTCTGGCCCTCCGTTCCCCACCTCTGCTCCA  
GCACCAGAGGACAGCACCCCTCTGCCAGGTGGACCAGGATGGAGGATCGCTTGACCAGGCGCCATTG  
CAGCCATCGTGATCGCTGCCCTGCTGGTACCTGCGTGGTCTGGCACTCGTGGTCGTCGCGCTAAGGAA  
GTTTTCCGCTTCT**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_028473
Insert Size:	366 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>Note:</b>	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
<b>RefSeq:</b>	<u><a href="#">NM_028473.1</a></u> , <u><a href="#">NP_082749.1</a></u>
<b>RefSeq Size:</b>	462 bp
<b>RefSeq ORF:</b>	366 bp
<b>Locus ID:</b>	73234
<b>UniProt ID:</b>	<u><a href="#">Q9CXL7</a></u>
<b>Cytogenetics:</b>	1 D
<b>Gene Summary:</b>	Plays a role in the regulation of chondrocyte maturation and postnatal endochondral ossification. May inhibit cell growth stimulation induced by FGF2.[UniProtKB/Swiss-Prot Function]