

Product datasheet for MC211404

Lrrc69 (NM_028499) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Lrrc69 (NM_028499) Mouse Untagged Clone
Tag: Tag Free
Symbol: Lrrc69
Synonyms: 1700034K16Rik
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC211404 representing NM_028499
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGCTGAAAGATTATTGGTAACAGCACTGAAGGGTGGCAAGAACAATAAGATCCTGACTTTGAATGGCA
 AGAGGATTACCAAGATGCCCTCGACACTGGAAAACTACCTAATCTGAAGACTCTAGACCTTCAAAAATAA
 TTCCATCTCCAAAGTATGCCAGAGCTAAGGACCTTGACCCAGTTGACGTTACTAAACCTGGGAAACAAC
 CATTTACAAGAAGTTCCTGAAGAGATAAAATATCTTACCTCCCTGAAGAACTTTCATTTATTTGGAACA
 GGATTTGCAGAATTGCTCCTGGAGTTTTAATGGCCTACATCGCTTGATCATGCTTAATCTGAATGACAA
 CAGACTTACGAGTCTTCCCTCAAGAAATTGGCAGGTTAAGGAGTCTTACCTACTTAAGCCTCAACCGTAAC
 AATCTCACCGTGATCCCAAAGAGCTTTGTTCCCTTGAGCATCTTCTGAACTGCATCTAAACTACAATC
 AGATCGGTACATCCCTGAGGAGATCAAGTTCCTGAAGAACCTTCAGCAGCTCTTCTGGTCAGGAACAA
 CATTGAAGAGTTGCCAGAGGAAATTTGTCATCTGAAAACTAAGAGTCTGGACATAGCTGGAATGTA
 ATCCAGATATTTCCGGCAGGCTTTCAGAACTGAGGCTGACAGAGTTCTACTGTGAGGGAACCCGCTCT
 TCCTGAAGAGGCCATTTTTGCTGTACAACCTAAAGACCTTGGACTCTACGGGAAATAGCAGCAAGATT
 TGTCTCTCTGTTAGAAAGAAAACGACCTTTGATTATGAATGTCATAGAGAGTTACCCAGAGGTCAAG
 GACAAACTCTCCAAAGCAAAAAAGTGTTCAATATGCAGGAAACCTTTTTAACTGAGTGGCTGGAGTGCG
 TTTACTTTGTGGCTCCATCGAAGAACTGGAAGATAAGCAGAAACCTAAAGCTGATTCTGTCCAGACATC
 AGTCTGTTCTATCAATGCTTTGACCAGCGTGATCCTGATGTCTTTGGAATTGCTCAAGAGTAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_028499



[View online »](#)

Insert Size:	1044 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_028499.2</u> , <u>NP_082775.1</u>
RefSeq Size:	1210 bp
RefSeq ORF:	1044 bp
Locus ID:	73314
UniProt ID:	<u>Q9D9Q0</u>
Cytogenetics:	4 A1