

Product datasheet for MC211419

Bpifa3 (NM_028528) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Bpifa3 (NM_028528) Mouse Untagged Clone

Tag: Tag Free

Symbol: Bpifa3

Synonyms: 1700058C13Rik; Splunc3

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC211419 representing NM_028528
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGCATCTACTCTGGAGGCTCCTGGTTCTCCTCCTCGGATTGCTAGCCTTGCCCTCAGCCCTGCCAAGC
 AGCCTTGGCCTGGACTGACCAAGGCTCACAAGGATGGCCGGTCCACCCTAGCGAGAATTATTGCACAGGG
 CCTCTGAAGCTCAATGCAGAAGGACGGATCCAGAGCATGCGCCTTCTGGACCGTCTGAACGTCTCAGGG
 ACAGTGGCTCCAGGGATGGTGGCTGGCTGATTGGTGGCATGAATTCCAGCAGCAACAAGAAATCAGCA
 TCAACATTACCAATGTGCAGCTGGACTGTGGTGGGATCCAGATGGCTTTCCCAAAGAGTGGTTCTCTGC
 AAACATCACACTTGAATTTGACATTGAATTCAAACTGCCCTTCAATTCGAATATCATCAAACGCATGCC
 TGCATGGGCCTCACTGCAGAGTCTGGCTGGAGAAAGATGAGTTTGGCCGGAGGGAGCTGGTGATGGGCA
 GGTGTGCAATGGAGCCCAGCAGTGGGGGTGCATCCATGTCCACTGAGGAAACCTCACCAGAGATGAAACA
 TTTTCTCCATAACCTCAGAGAAAGCCTAGGGAAGTTATCCCAAATCTGGTAGAAAGTCAGGTCTGTCTC
 CTGATTGGTGAGATCTCCGGCAGCTGGATGTGAAGCTGCTGAAAGGCCTAGTAGAACAGGTGCCTGCTC
 ATAATCTGGGTCAACCC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_028528

Insert Size: 720 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_028528.3</u> , <u>NP_082804.1</u>
RefSeq Size:	1110 bp
RefSeq ORF:	720 bp
Locus ID:	73388
UniProt ID:	<u>Q9D9J8</u>
Cytogenetics:	2 H1