

Product datasheet for MC211429

1700074P13Rik (NM_028550) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 1700074P13Rik (NM_028550) Mouse Untagged Clone
Tag: Tag Free
Symbol: 1700074P13Rik
Synonyms: AV047902; Tryx5
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC211429 representing NM_028550
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAGACTTTTTTCATCTTCACTCTCTTAAGTCTCGCAGTAGCTAATACTCCCAAAGCCGTTCTGAAAG
AATATAACAATTCAGAAGAATACTTACCAGAGAATTTTAAATGTTCTTATATGGTCTATCTGCAGTCCAG
CCCAGAGCCTTGTGTTGGGACTCTCATTGACCCTCAGTGGGTGCTCACCGCTGCTCATTGTTCTTTACCA
ACTAAAATCCGACTGGGAGTTTACCGACCAACATCAAAAATGAGAAAGAACAGATATGTAATTACTCAT
TCACTGTAGTCCATCCTAATTTTATGATGCAAAATTAAGTGAATAATGACTTGATGCTGATAAACTCTCATA
TCCTGCTACTATCAACATGTATGTGGAACTATAGCCATAGCTATGGAACCAATGGCATTAAATGAATCC
TGCTTCATACCAACATGGACCTGGAATAACTACAAAACCTTAGTGATCCTGACATCTTGACATGGATAA
ATGAGTATTCTTTCCCAAGTGACTGTTTGGACTCTTCAACAGAAACAAGAAACAAGAATAAA
CATTATGTGCATAGGACATTCACCTAATGCTATGTCTGCAACTAAGGAAGTATCAGCTGCTCCAGCCATC
TGCAGTGGGAGGGTGCATGGAATATTGCTCTGGGGAAAGCAAGTGTAGCCAATGGAAGCAAAGGCTTCT
TCACTGAAATTCATCCCTATGCAAGATGGATCTTAAAAATGATACGATCACACTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_028550
Insert Size: 756 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_028550.3, <u>NP_082826.1</u></u>
RefSeq Size:	1647 bp
RefSeq ORF:	756 bp
Locus ID:	73481
Cytogenetics:	6 B1