

Product datasheet for MC211457

Lce1l (NM_028628) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Lce1l (NM_028628) Mouse Untagged Clone
Tag: Tag Free
Symbol: Lce1l
Synonyms: 1110008K04Rik; AI415426
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC211457 representing NM_028628
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCCTGCCAGCAGAGCCAGCAACAGTGCCAGCCTCCTCCCAAGTGCACCCCAAGTGCCTCCCAAAT
GCCAGACTCCTAAATGCCCCCAAAATGTCCCCTAAGTGCCCTCCCAAGTGCCTCCTGTCTTCTCTG
CTGTAGCCTGGGTTCTGGGGGCTGCTGTGGCTCCAGCTCTGGTGGTTGCTGCAGCTCTGGTGGCTGCTGC
AGCTCTGGTGGTTGCTGCAGTTCGGGGGTGGAGGCTGCTGCCTGAGCCACCACAGGCCTCACAGATCTC
TTCGTCGCCATCGACACAGCTCTGGATGCTGTAGCAGTGGTGGTAGCAGTGACTGTTGTGGCAGCAGTGG
TTGCTGTGGCAGCAGTGGCTGCTGTGGTGGTAACCAGCAGTCTGGGGACTGCTGC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_028628

Insert Size: 408 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_028628.2](#), [NP_082904.1](#)

RefSeq Size: 670 bp

RefSeq ORF: 408 bp

Locus ID: 73730

Cytogenetics: 3 F1