

## Product datasheet for MC211482

### Them7 (NM\_028747) Mouse Untagged Clone

#### Product data:

Product Type: Expression Plasmids

Product Name: Them7 (NM\_028747) Mouse Untagged Clone

Tag: Tag Free

Symbol: Them7

Synonyms: 0610012H03Rik

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC211482 representing NM\_028747  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGC**

ATGGCCATGAGGCGGCTGCTGAGAAATATCTGTTTACCTTGAATCACAGGAATTCTGACGCAGCAATCA  
CTATGATTCTTGGGCATGGAGGCTTCTCCTAGAAATGCCGCTCGTTCGTATCTACATCATCAGCCAA  
TCAGCGAGCCGAAGATCCCCACCGGCAGCTGGAGGCCTATGGATACTTCTCCCATCCAGACCCGGTGG  
CAGGACAATGACCAGAATAGCCACGTGAACAACGCTGTGTACTACAGCTACTTCGACACTGTAATCAGCC  
ACTACCTGATCAGGTACTGTGGCCTGAAAACAGCTTGCTGACATCCCCCAGGTGGGGTTCATGGTGAC  
CAACCAGTGCAAGTTCCATACACCCATCAGCTTCCTCAGACCCCACTGGCTGCCCTGGCTGTGGAGAAA  
GTCGGTCGCTCCAGTGTGTGTACCGACTGGCTCTGTTTCCTCCTAAACCCACCAAGGAGCTGCCTTCAG  
TCAATCACCATGATCTCATTGATGGTTTGTTCGAGGCCACCCTAAATTAGCACAGTTTGACACTCTGGC  
ATGTACCACGGGGTCTTCAGTCCATGTCTTTGTGAATCCAGCCACTAACAAACCAGAGAACCTTCCAGAG  
GATTTCAGGAGGGGTCTCCAGAAGCTGATGAGTCCAGCA**TAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM\_028747

Insert Size: 672 bp



[View online »](#)

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).                                                                                                                                                                    |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                              |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>Note:</b>                  | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>RefSeq:</b>                | <u><a href="#">NM_028747.2</a></u> , <u><a href="#">NP_083023.1</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>RefSeq Size:</b>           | 1598 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>RefSeq ORF:</b>            | 672 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Locus ID:</b>              | 74088                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Cytogenetics:</b>          | 2 E3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |