

## Product datasheet for **MC211492**

### **Fbxl22 (NM\_175206) Mouse Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Fbxl22 (NM\_175206) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Fbxl22  
**Synonyms:** 1110004B15Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >MC211492 representing NM\_175206  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGCACATAACCCAGCTCAATCGGGAGTGTCTGCTGTGTCTCTTCTCCTTCTGGACAAGGACAGCCGGA  
GGAGCCTCTCCAGGACTTGCTCCAGCTCCGAGATGTGTTTGAGGACCCACCCTCTGGCCCTGCTGCA  
CTTCCATTCTCTCGCAGAGCTCAAGAAAGACAACCTCCGGCTGAGCCCTGCGCTGCGCAGCCTGTCCATC  
TGTTGGCATTCCAGCCGTGTGCAGGTGTGCAGCATCGAGGACTGGCTCAAGAGCGCCCTCCAGAGGAGCA  
TCTGCAGCCAGCAGAGAGCCTGGTTAATGATTCCTCCTGCAGGTGTGCAACAGGTGCCCAACCTGAC  
GTCCGTACGCTGTCCGGCTGCGGCCAGGTACGGACGACTGCCTGGCGCGCTTCTGCTCAGCTGCCCG  
CGCCTGCGTACGCTGCGCCTCGAGAAGTGCAGCGCGCTTACCAACCGCACGCTGGCGGCCGTGGCCGCGC  
ACGGGGCGCGCTGCAGACTCTGCACGTGGACTTCTGCCGCAACGTGAGCGCCGCGGCCCTGCTCCGCT  
CCGCGCCGCTGCCCAACCTGCGCCTCAGCGCCGAACGAGCGCGGCCATGATACCCGATCAGCCGCCG  
CGCGCACGCGTGTGCGCTGCCTCAGCTTCTCGGCTGCACCGGGAGACAAAGCATGCCTACCCACGG  
GCTCCTACTGA

AGCGGACCGACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
TGGATTACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-RsrII  
**ACCN:** NM\_175206  
**Insert Size:** 711 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_175206.4</a> , <a href="#">NP_780415.1</a>
<b>RefSeq Size:</b>	3935 bp
<b>RefSeq ORF:</b>	711 bp
<b>Locus ID:</b>	74165
<b>UniProt ID:</b>	<a href="#">Q8C7B6</a>
<b>Cytogenetics:</b>	9 C
<b>Gene Summary:</b>	Substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex. Promotes ubiquitination of sarcomeric proteins alpha-actinin-2 (ACTN2) and filamin-C (FLNC) (By similarity).[UniProtKB/Swiss-Prot Function]