

Product datasheet for MC211529

Arl13a (NM_028947) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids

Tag: Tag Free

Symbol: Arl13a

Synonyms: 4933422M2IRik; Arl13

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC211529 representing NM_028947
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTTCCGACTTTTGACTGCCTGCTGGTCTCGTCAAAAGACCACAGAAGAAAAACAAGGAATGTGACTA
TCATTGTTATTGGTTTGACAACCTCTGGCAAAAGTCGTCTCATAGAGGCCTTCCAAAGACTAATTCCTAG
CAAGATGCACAAAGAAATGAGACCAACACAAACACGCTTCTACTAGATGATTACCAAGTATCTATCTAT
GACCTAACTGGAGACGTAAGGCGCTGAGAAATGGTCAAGTTACTATGCTCAGGCCCATGGGCTTATTT
TCGTGGTGGATTCTAGTGACATAGCACGCATCCAGGAAGTGAAGATCATCTTAACCTCGCATAATGTTTGA
CAAAAGGGTCTCAGGGAACCCATCTTAATCCTGGCAAAACAAGACAAGAAGAACGCTCTTTTGCCT
TGTGATATCATTGAATATCTACTTCTGAAAGGCTCGTGAATGAGACGAGATCCATGTGCCGCTGGAGC
CATGTTCAACCATAAACGCATTCCAAAGAGGCAGCAAAAGCCAATAGTTGTAGGGCTGCGATGGCTACT
AGCTGCCATTGGGGATCAATACAGTGAGCTGTATACTCGACAACAAATACCCAATATGAGCATCTCAACC
TCCAAGAACATCAGAGGCTGTGGAGAAAGATGCTTATCAGACAGCTTACCTAACAGAGTCAGACACAGCA
ACAGACAGCATTTTGATAAAAGACAACACATGGAGCATAGACAGCATGTGGAGCAAAGGCATTTTGAGAA
AAGACAGCATTTTGAGGCAAGACAGCGCTTAATACAACGTTCTCTGGATGCTAGGCCCTAAAGCCAATT
CTACAGAAAGATGGTTTTCGAATAAGGCCCAAAAAGAATATGTCAGTAACATTTGCTCTAGATGTCAATCA
TGGAGGAAGGTGAATGTTCTCGGAAGATTGGAGCTCCAACCATCAGTAAGCCTTGCAAGTGTAGTGTTA
TGATACAAAGACTCCAGACCCATCTGTTGATGTCAGTCTTTTCAAAGCTCATAAGCCCAAAAGAGAATA
GAAACCTGGGACACTGAAGAAATGTTCTTGAAGACCTCTGTAAGAAACCTTCCATTCTGTGCAATAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



ACCN:	NM_028947
Insert Size:	1119 bp
OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery. The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_028947.1 , NP_083223.1
RefSeq Size:	1742 bp
RefSeq ORF:	1119 bp
Locus ID:	74448
UniProt ID:	Q9D416
Cytogenetics:	X E3