

Product datasheet for MC211536

Elov17 (NM_029001) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Elov17 (NM_029001) Mouse Untagged Clone
 Tag: Tag Free
 Symbol: Elov17
 Synonyms: 9130013K24Rik; AI840082
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 Cell Selection: Neomycin
 Fully Sequenced ORF: >MC211536 representing NM_029001
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGCGTTCAGTGATCTTACATCGAGGACTGTGCGTTTTTATGATAATTGGATCAAAGATGCTGATCCGA
 GAGTTGAAGACTATCTCCTCATGTCCTCGCCTCTGCCACAAACCATCATCTGGGCCTCTATGTCTATTT
 TGTTACGTCCTGGGACAAAGCTCATGGAGAACCGGAAGCCCTTGAAGTCAAGAAAGCGATGATAACG
 TACAATTTTTTTCATAGTACTGTTTTCTGTGTATATGTGTTATGAGTTTGATGTCTGGCTGGGGTACAG
 GTTACTCGTTTTGATGTGACATTGTTGACTATTCTCAGTCGCCAAGAGCAATGAGGATGGTGACACCTG
 CTGGCTTTATTACTTCTCCAAATTCATCGAGCTGTTAGACACTATCTTTTTTGTCTGCGTAAGAAAAAC
 AGCCAAGTGACTTTCCTTCATGTCTTCCATCATACGATCATGCCGTGGACCTGGTGGTTTGGAGTCAAAT
 TTGCTGCAGGTGGTCTGGGAACATTCCATGCCTTTTTAAATACAGCTGTGCACGTGGTCATGATTTCTTA
 CTATGGTTATGTGCAATGGGACCAGCCTACCAGAAGTATTTATGGTGGAAAAACATTTGACATCTTTG
 CAGCTTGTCCAGTTTGTCTTGTCACTATCCACATAGCCAGATCTTCTCATGGAGGACTGCAATTACC
 AGTACCCAGTTTTTCTGTATATCATCATGAGTTACGGATGCATCTCTGCTGCTTTTCCACTTTTTG
 GTACCGTGCTTACACCAAGGGTCAGAGGTTGCCAAAACCTAGAAAATGGAAATGCAAAGCAAGCGC
 CACTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
 ACCN: NM_029001
 Insert Size: 846 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_029001.5 , NP_083277.3
RefSeq Size:	5453 bp
RefSeq ORF:	846 bp
Locus ID:	74559
UniProt ID:	Q9D2Y9
Cytogenetics:	13 D2.1
Gene Summary:	Catalyzes the first and rate-limiting reaction of the four reactions that constitute the long-chain fatty acids elongation cycle. This endoplasmic reticulum-bound enzymatic process allows the addition of 2 carbons to the chain of long- and very long-chain fatty acids (VLCFAs) per cycle. Condensing enzyme with higher activity toward C18 acyl-CoAs, especially C18:3(n-3) acyl-CoAs and C18:3(n-6)-CoAs. Also active toward C20:4-, C18:0-, C18:1-, C18:2- and C16:0-CoAs, and weakly toward C20:0-CoA. Little or no activity toward C22:0-, C24:0-, or C26:0-CoAs. May participate in the production of saturated and polyunsaturated VLCFAs of different chain lengths that are involved in multiple biological processes as precursors of membrane lipids and lipid mediators.[UniProtKB/Swiss-Prot Function]