

Product datasheet for MC211551

Lrrc74b (NM_001145435) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	Lrrc74b
Synonyms:	4930451C15Rik; mCG12354; mCPI7390
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC211551 representing NM_001145435 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGAAGGGTCCCTGCGAAGTTGAGAAGAATGAAGATCAAGAGGGAGAAGCTGCGGCCACTGGGCCTCAGG
CAGAGACCCCTGGAGGCCGAAAGGTCCTGGACTGCAGATTACACTCTGCTCTGGAAGCAGAAGGCACCCA
TGGGCTTGGAGAACGGGTATGGCTACCTCTATCTGAAGTCTTGCCGGGCCAATAGTGTGTGCCAGTT
TCTTGTCTTCTGCGCCAAGAGGGTGCCTCTGAATTGAACCTTCGGCACCGTGGCCTGGGGCCCCAGGGGG
TCCGGGCATTGGCTTCCGTATTGACCTCTAATCCCTACATCAAGCGTCTGGATCTTCGGGACAATGGGCT
CTGTGGTGCTGGAGCAGAGGCTCTGGCTGATGTCTGCGTAAGAACAGCATTATCTCTGATGTGGACCTG
TCAGAGAATCAGATTGGAGCTGCAGGGCTCCAGGCGATCTGCACTGCCCTTGCCTTAAACCCAACGTGTG
AGAAGATGCAGTGCAGGGAACAGGCTGGAGGAACAAGCAGCCCAGCACCTCGCTGCTCTTCTGCTGCA
TCACAGGGGCTGAAATCTCTAGACCTCAGCTATAATCAGCTCAATGACCTGGCAGGAGAGATACTTGA
CCTGCTGTGGCAGAAAACAGGACTCACAGAGCTTAACCTCAGTTGGAATCACCTGCGAGGTCTGGGAG
CCACCGCGTTTCCCGAGGCTGGAGGCAAATATTTCTTGAAGTCTAGACATCTCTCATAATGGCTT
TGGAGACTCTGGTGCCAGTGCGATTGGTGATGCCCTCAGGGTCAACAATGTTTGAAGAAGTGAATATG
AGGAACAACCGAATCTCCGTGTGAGGAGCCCTGAAGTTGGGATTGGGTCTGCAAGTCAATCAGACGCTGC
GGATTCTTATTGATATCAAGTAAGCAGAGAATGTGAGGACCTCGTAGCTCCATGAGTGAATTTCTTCC
AGGACTCTGTATAAAGCGTTATACTTCTAGAAGAAAAGACTGGCCACAAGCCTCTACACCATCCAGCCA
GCATCTGCACCAAGTGACTCTGGTCTCAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



ACCN:	NM_001145435
Insert Size:	1080 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001145435.1</u> , <u>NP_001138907.1</u>
RefSeq Size:	2135 bp
RefSeq ORF:	1080 bp
Locus ID:	74685
UniProt ID:	<u>Q14BP6</u>
Cytogenetics:	16 A3