

Product datasheet for MC211554

Pih1d3 (NM_029062) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pih1d3 (NM_029062) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Pih1d3
Synonyms:	4930521A18Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC211554 representing NM_029062 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGAATCAGAAAATGCTGAAGCTGAGAATGTGGTAACAGGGAATTTGGAAGCAAAGAGCAAGGAAGGTA
AAAATATGCAATTTGTGACATTGTCTTCTGCCTGGTCCCTTCAAGCACTATCCAGTCTCCTCAATCCCGA
AGAAGAGGATGGTTTTGATTTTGAACAGGGGCAATGCTCATTTACTATTGGAGCCATGGGCCCTGGGAAT
ATTGGACCACCTAAAGCAAAAGAGTCAAAGCTATTCCTGAACCCAGGAGTGATGAAAGTGAGAACATCT
GGAATCCAGAAGAGGTTCCAGAAGGAGCAGAGCATGATGATATATGGGATGTTTCGTGAAATTCCAGATTA
TGAGATTGTATTCCAACAACTGTGGAACTGAGGACGTATACCTAGGATTAACAAGAAAGGATCCTTCA
ACGGCATGCTGCCAGGAGCTGGTGGTGAATAAACTGCCCAATACAAACCTTCTGAAATTCAAATTG
ATATCCAAGAAATGCTCCTTGATCTTCGCACTCCTAGAAAGAAGCTGTTAGTAACTTTCCCCAGCCCGT
GGAACGGAACAGTGCTAGAGCATCCTATATCTGGAAGCTGAAACACTTGAAGTCAGGATGACCGTCCAA
AGGGATTTGGATTTTAATATTTCC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_029062
Insert Size:	657 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_029062.2 , NP_083338.1
RefSeq Size:	1449 bp
RefSeq ORF:	657 bp
Locus ID:	74708
UniProt ID:	Q3KNI6
Cytogenetics:	1 B
Gene Summary:	Plays a role in cytoplasmic pre-assembly of axonemal dynein.[UniProtKB/Swiss-Prot Function]