

Product datasheet for **MC211573**

4930503E14Rik (NM_029131) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 4930503E14Rik (NM_029131) Mouse Untagged Clone
Tag: Tag Free
Symbol: 4930503E14Rik
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC211573 representing NM_029131
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTTTTCTGGCTGCGCAGGCTATTTACAGGGAGAATGGCGATCAAGGAGAGACAAGACCAAGGCAGA
AGGAATCTGCCATCCCTTCTGTAAAAACAGAAGAATGAAATCATTCTGGAAAGGCACATGTCTGCTGG
GAAAACATCATCCAAAATTGCAACATCACTAATCATATGAAGAATATGAATAAACTAGATGATGAAA
TTTTATATCAGAAAGATAAATGCTGAGAGGTTAGAAGCTTTTTCGAATCTGGACATTGACATGAATACTG
ATTTGAACTACAGGATGAACATTGAATTCACCATCATCAATCACAACATGAGAAGACAATGTTGGATAT
GGAAAAATGACGCAGTCAATAAGTGATACCATTGAGAAGTACAAGGAATTCATAGAAGATAAGGATTCC
TACAGGTGTGAAATACTCCAGCAGAACTGGAACATGGAACAGACCAGGACAAGGTCTCCTTCGAAGAGA
ACATCTGCGTGTGCGAAGCAAATATTTTCATCTCATCCCCAGTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_029131
Insert Size: 534 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_029131.3](#), [NP_083407.1](#)

RefSeq Size: 1289 bp

RefSeq ORF: 534 bp

Locus ID: 74954

Cytogenetics: 14 C1