

Product datasheet for MC211596

Tmem269 (NM_029198) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tmem269 (NM_029198) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Tmem269
Synonyms:	4930538K18Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC211596 representing NM_029198 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCTTGGCACAGGAGACCCCAATCTGCACTCAAGTTGCATCCTCTCTAGGAAAAAGGAAGACAAGGTC
 ATCACATCATGAAACATGAGCAGATCGGTTCTCAGATGACTGAATTTCTGGAAAGACATTATTAATGC
 CTTGGCCCTGGCCAACATGATCCTGGGGCTCTTCTCCATCTTCTGCAGCTTCAGCAGGAAATCCTACTGT
 GCTTCCTGGATGCTTCTGATCAGCTTCTGTTAGACATTGCAATAGGAACGATGACCAAAACCTCAACA
 TCCCACAAACTGGGACTCGAAGTGAATGATTTTGCATCTTACCACCTTTGGCCTGGCCTCGGCTCT
 GCTGCTGGGCGTGGATGGGCCTCTGAATGGGTTCTGGCCATCATCTATGTGCTAACCACCTCTTTCCGC
 ATGTGTTTTTATTCAACTGGAGGAGCTACCTCTGGATACAAGGGCCTACCTGCCCTATGCTTCCTGTG
 TTTTGGCCTCCACCTGCCTTCTGACTAAGGGCAATACATTCATTCTCTGCTGCATGGCTTCACTCATGAT
 TCTCTTCATGATCGACCAGAGTTGCTACCCACACGATGAAATCCTAGATTCAGACAACTGGAAAAAATC
 GTTTATATAGGAGGTGTATCCTGCTGTTTTCTCTCCGTTCCCGCTCACCAGCTTTTACTGCCTGACCT
 GGTGCTCTCTACATCTTCTCCCTGAAACCTGTGGGCAGAGGAGTTCGTATTAAGCCCTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_029198
Insert Size:	765 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_029198.3</u> , <u>NP_083474.1</u>
RefSeq Size:	1504 bp
RefSeq ORF:	765 bp
Locus ID:	75180
UniProt ID:	<u>Q9D4Y8</u>
Cytogenetics:	4 D2.1