

Product datasheet for MC211601

Spaca6 (NM_001162909) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Spaca6 (NM_001162909) Mouse Untagged Clone
Tag: Tag Free
Symbol: Spaca6
Synonyms: 4930546H06Rik; 9530079D04Rik; B230206P06Rik; let-7e; mir-99b; mir-125a; Ncrna00085
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC211601 representing NM_001162909
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGACCTCACAAGGTCCTACTATCTTCTCCTCAGACCCGGAGACCCTCTGTTATGGGTCTTATTTCTTAG
 TAGGGAGCATTGCTCTGTTGTTCTCTGCTGATCTTCAGGGCATCCACGTGGGCCTGTCTCTCTGCTTAC
 TACCTATGAGGAACGCCCTCCGTGTCTGCCAGCTGTTTGTGGGCAGAGAGGAAACCAAGATAAACCTGTGC
 AGAAATGAGCTTGAAAGGGCTTTTGAAGATCTCAAAGACATGAAAATCAACTATGATGAGAGGAGCTACC
 TGCATGATGAATTTACACAGATGACAGTTTCTCTTCAGGAGAAGGCTGCTAGACGACGGGAGCCCTTTTG
 GCTTGCCTTCAAAGATGCTGCTGCAAAATTGAAGAGGACCATTGAACACCTTAAAAAGCCCGGCTGC
 ATCCCTCCTTGTGGTCTCCAGGAGTGGCTAGGCTTTTCCACTGCTCGGGCTGCTTTTCTAAACTCTGCG
 ATCTGCCTCTCGATTGTCCAGTTCAGGACATGCTGGTGAACCGAGGGGACCAGGCTTTGTTTTCTTGAT
 CGTGGCTTTCGAGTTGCCAGAGTCGGAGATCACTTATTCCTGAAATTTGTGGGAGGTGTCCGAACCTAAG
 GACGTTACGTAATCCGTGATATGCCGGGGCTCACGGTACCTGGCACGAATACGGCCAGTGAACCCCA
 AACACGGTGGGACCTTTTCTGCGTGATCCTGCACGACGACGCCCCCTGGCGCGGCTCTACTTTTATCT
 GAACGTGACGGGCCCTCCTCCGCTGAGGACACTGAGCTCCAGGTACGTTCCGGGAAGTATGAATCGG
 ACGCCGCGGAGCCGAAATGATCCAGCCCTGGAGCCCCAGCCTAGGCGAACTACTTACCAACCCCAAGG
 CTCTGACGCTTGCAACCTATTCCTGCTCGCGGCCACCGCTGCACTTGGGTCTGCCAGTGTACTCTGCT
 AGTGTGGCTGTTTTTTCGATGTTACTTAAGTGGCAACT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001162909



[View online »](#)

Insert Size:	1020 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001162909.1 , NP_001156381.1
RefSeq Size:	1102 bp
RefSeq ORF:	1020 bp
Locus ID:	75202
UniProt ID:	E9Q8Q8
Cytogenetics:	17 A3.2
Gene Summary:	Sperm protein potentially involved sperm-egg fusion.[UniProtKB/Swiss-Prot Function]