

Product datasheet for MC211631

Cklf (NM_029295) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cklf (NM_029295) Mouse Untagged Clone
Tag: Tag Free
Symbol: Cklf
Synonyms: 1700001C14Rik; 1810018M11Rik; A730028G07Rik; C32; CKLF1; Cklf2; CKLF3; CKLF4; CKLF5
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC211631 representing NM_029295
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGAGACTCCACGGCCGGTCGTAAGCCGTCGGCCCTTCTGCTGTACTCTGAAATGCTTTGTGAAGTTCC
 TCGGGTTGGTGGTAACTGTGACATCTATGATCTTTTTTCATCGTGGGACAAGCCCTGAACCATATATCGT
 CATCACTGGGTTTGAAGTCACCGTTATTTTCTGTTTCTAGTCTTATACACGTGTGGACTTGACAAGATA
 ATGAGATCTTTCTTTGGCCTTTGCTTGATGTTATTAACCAATGGTAACAGCATTATGCATGCTCATTG
 TGCTGTGTTGGCTCTAATACCAGAAACCTCAACAAAGACAATCCTTGGAGGGGTGTTTGGTTTCTGAC
 AGTAACATGTAATAATGCCGACTGTGCCCTTATGTGCCAGAACTTCGGTTTCGTCCACGTCAACCTTAT
 CAGAAAAAGTCTACTAATGACATCGATGACAGAGAA**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: NM_029295
Insert Size: 459 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_029295.3</u> , <u>NP_083571.1</u>
RefSeq Size:	667 bp
RefSeq ORF:	459 bp
Locus ID:	75458
UniProt ID:	<u>Q9DAS1</u>
Cytogenetics:	8 D3
Gene Summary:	May play an important role in inflammation and regeneration of skeletal muscle. [UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest isoform (1).