

Product datasheet for MC211637

Spata19 (NM_029299) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Spata19 (NM_029299) Mouse Untagged Clone
Tag: Tag Free
Symbol: Spata19
Synonyms: 1700001L23Rik; AV047838; spergen-1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC211637 representing NM_029299
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGATCATTACAACATGGATTATGTACATCTTTGCCCGAAAACCGTAGGGCTACCCTTCCCACCAAGGG
TGAACCTCGGACATTGAAGTTGAAGAAAAGCGAAGCTGTATCCGTTGTGCAGCACTGGTTGAACAAAACGGA
AGAGGAGGCTTCTCGGAGCATAAGGGAGAAGATGTCCATCAACGACTCTCCACCCATGGACATGACATA
CATGTGACCAGAGATTTGGTGAAACACCATCTCTCAAAGTCTGATATGTTAACAGACCCGAGTCAAGAAG
TCCTGGAGGAAAGAACAAGGATCCAGTTTATAAGATGGAGCCACACCCGTATCTTCCAAGTGCCAAGTGA
AGTGATGGATGACGTCATGCAGGAACGAATAGATCAAGTGAGACGAAGCGTATCTCATCTCATGTGTGAC
TCGTACAATGATCCAAGCTTCAGGACTTCTTGCTCAGAATGCTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_029299

Insert Size: 465 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_029299.3</u> , <u>NP_083575.1</u>
RefSeq Size:	711 bp
RefSeq ORF:	465 bp
Locus ID:	75469
UniProt ID:	<u>Q9DAQ9</u>
Cytogenetics:	9 A4
Gene Summary:	May have a role in spermiogenesis.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).