

Product datasheet for MC211655

Paqr9 (NM_198414) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Paqr9 (NM_198414) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Paqr9
Synonyms:	1700020G04Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC211655 representing NM_198414 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGCCGCGGCGCCTGCAGCAGCGGGGCGCGCGTGAAGGGACCGCGGCCTCGACCTCCCGGCGCTCGC
ATCCCGCCTCCGCTCCGCCCCCGAGCCCGCGCGGCCACGACCAAGCCGCTGCTGCGCTGGGACGA
GGTGCCCGATGACTTCGTGGAGTGCTTCATCCTGTCCGGTACCGCGGCTGCCGTGCACCGCGCAAGAG
TGCTGGCCTCGGTGCTGAAGCCACCAATGAGACGCTCAACTTCTGGACGCACTTCATCCCGCTGCTGC
TGTTCTCAGTAAGTTCTGCCGCTCTTCTTCTGGGCGGCAGCGATGTGCCCTTCCACCACCCCTGGCT
GCTGCCGCTGTGGTGCTATGCGTCGGGGGTGCTGCTGACCTTCGCCATGAGCTGCACGGCGACGTGTTT
AGCTGCCTGTGCTGCGCTGCGCGCCGCTTCTTCTACCTGGACTACGCTTCCATCAGCTACTATGGCT
TCGGCAGCACCGTGGCTACTACTACTACCTGTGCCCAGCCTGAGCCTGCTGGACGCCAGGGTCATGAC
GCCGTATGTGCAGCAGCGCTGGGCTGGCACGTGGACTGCACCCGCTCATCGCGGTCTACCGCGCGCTC
GTGCTGCCCGTGGCCTTCGTGCTGGCTGTGGCCTGCACCGTGGCTTGTGTAAGAGCCGACCGACTGGT
GCTCCTACCGTTCCGCTGCGCACCTTCGTCTTCGTATGCCGCTCAGCATGGCCTGCCCATTTATGCT
AGAGAGCTGGCTCTTCGACCTGCGAGGCGAGAACCCACACTCTTCGTGCACTTCTACCGCGCTACTTC
TGGCTGGTGGTGGCCGCTTCTCAACGTGAGCAAGATCCCTGAGCGCATCCAGCCGGGTCTCTTCGACA
TCATTGGCCACAGCCATCAGCTCTCCACATCTTCACTTTCCTCAGCATCTATGACCAGGTGTACTACGT
GGAAGAGGGCCTGCGCCAGTTCTCCAGGCACCCCTGCCGCGCCACCTTCTCGGGCACAGTGGGCTAC
ATGCTGCTGCTGGTGGTTTGCCTGGGACTAGTCATCAGGAAGTTCCTGAACAGCACCGAATTCTGCAGTA
AAAAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Chromatograms:	https://cdn.origene.com/chromatograms/ja1866_g11.zip
Restriction Sites:	Sgfl-Mlul
ACCN:	NM_198414
Insert Size:	1128 bp
OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery. The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_198414.2 , NP_940806.2
RefSeq Size:	2305 bp
RefSeq ORF:	1128 bp
Locus ID:	75552
UniProt ID:	Q6TCG2
Cytogenetics:	9 E3.3
Gene Summary:	Plasma membrane progesterone (P4) receptor coupled to G proteins. Seems to act through a G(s) mediated pathway. May be involved in regulating rapid P4 signaling in the nervous system. Also binds dehydroepiandrosterone (DHEA), pregnanolone, pregnenolone and allopregnanolone.[UniProtKB/Swiss-Prot Function]