

Product datasheet for MC211686

Idnk (NM_198004) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Idnk (NM_198004) Mouse Untagged Clone
Tag: Tag Free
Symbol: Idnk
Synonyms: 5133401N09Rik
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC211686 representing NM_198004
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGAGGCACCCGGGTGCTGCTGGTGATGGGTGTGAGCGGCTCCGAAAATCCACTGTGGGTGCGCTGC
TGGCCTCCAAGCTGGGATGGAAATCTATGATGCCGATGATTACCACTCCGAGGAGAATCGGATAAAGAT
GGCGAAAGGGGTACCACTGAGCGACCAGGACAGGATCCATGGCTTTGCACCTTGACGACATTTTACTA
AGAGATGTGGCTTTGGGACAGCCTGTTGTTCTAGCCTGTTGAGCTCTGAAGAAAACGTACAGAGACATCT
TGATCCGAGGAGGAAGTATGCGCCCTGAAAAGTATGACTCAGCAAAGGAACCGCTGGCTGGTGGAA
GCTCCTGGTGGTCTACCTCTGCGGCTCGTTTGACATCATTTATGGACGCTTGCTCCAAAGAAAAGGACAT
TTTATGCCGCCGAGTTACTGCAGTCCCAGTTTAGTATTCTGGAGCCCCATCAGCTCCCGAAAACCTCA
TCCAAGTCAGCGTGGACAAAAGTCTCCCGGAGATCACCCTGCCGTTATGGAAGCCCTCAAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_198004
Insert Size: 555 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_198004.3</u> , <u>NP_932121.2</u>
RefSeq Size:	1511 bp
RefSeq ORF:	555 bp
Locus ID:	75731
UniProt ID:	<u>Q8R0J8</u>
Cytogenetics:	13 B1