

Product datasheet for MC211733

Fam187b (NM_175240) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fam187b (NM_175240) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Fam187b
Synonyms:	1700020B09Rik; Tmem162
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC211733 representing NM_175240 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGCTGGCTACTCTGTGGCTGGTTGGTCTTTCTCTTCCCATGCTGTGGGCCAGAGACTAATCAGCTGTC
 CCTACAAGAATGTGTGTCAGTATGCCCTGCTCTCAGGCAGTGACGTGATACTGCACTGTAACACCCAA
 AGCACTTTGGTATTTTCTTCATCTTTAGAGGACAACTGTCCCTGGTCAACTCCAAGCCTGATGGGAGG
 GTACTGCCTGGGAGCGACCTGCAGTTGTCTGACCTAAGCCTTCCAGACAGGCCTCTACCGCTGCCTGG
 ATAATCACAAGGCCCGTTTGGTAGAATATGAGATTGACTTTCAGAACATAGCTCTGTTGCACATCACTCA
 CAAAGACCTGGGCCAAGAGCCCATGGGAAATGAGAGCATGAACCTGGGAGGCAAGGTTCTTGTTTTTACC
 CGCTGGGATCCTTGGCAAGATTGTAACCGTTGCCAAAAGCCAGGTGAGCGCAAACGCTCTAGGCTACTGCT
 ATGTGGAGGAGCCCAAGAAAAGCCCATGCCTTGCTGGCTCTACCTGAGAGAGGAGAAGGTGACCAACAG
 CCGCTGAGGCCTGAACCTCCAGCTACAAGCCTGCCAGGTACCTGTGACACCGCCACGGAGACCAACAG
 CCGTACTTTGTCTTTGACACATATCTGCTGGACAAACCAACAGCAATGCTCGACTCAAGTGTCCCTTGG
 CTTCCATCTACAGGCCTGTCCACTGGGAAGCTGACAACAGCCCTCTGACATGGCAGGACCAAGCTCTCTGG
 CCAGACGGTCAGCACCATCATGGACCTGCACAGTGGTGGCCAGCACCTGAAGGTCTTCCAGCCTGCCACT
 TACAGATGTTTTGTGGAGCAGGAGCTCATAGCCAGTTCAACCCCACTGAGAGACAATCCAAAGCCGAGA
 ACCCTTGGCAGCCAGAAATCCAGCCAGACAAGGCAGACTCCGTTCTTAGGAGATTGAAGTTGATGGTGCT
 GAGTATCTCTGTACTGGCCGTTGGAGGGCTTCTCTGTAAGGTGGTTTTTCGTCCTGTCTGCGCAAGAAG
 AGGAGTCAGGTCCTCTGGTGAA**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_175240
Insert Size:	1077 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_175240.4</u> , <u>NP_780449.2</u>
RefSeq Size:	1485 bp
RefSeq ORF:	1077 bp
Locus ID:	76415
UniProt ID:	<u>Q0VAY3</u>
Cytogenetics:	7 B1