

Product datasheet for MC211739

Krtap4-7 (NM_029613) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Krtap4-7 (NM_029613) Mouse Untagged Clone
Tag: Tag Free
Symbol: Krtap4-7
Synonyms: 2310037K05Rik; KAP4.7; KRTAP4.7
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC211739 representing NM_029613
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGTCAGCTCCTGTTGTGGCTCTGTCTGCTCTGAGGAGGGCTGTAGCCAAGGCTGCTGCCAGCCAGCT
 GCTGTGTGTCCAGCTGCTGCAGGCCTCAGTGTGCCAGTCTGTGTGCTGCCAGCCACCTGCTGCCGCC
 CAGCTGCTGCATCTCCAGCTGCTGCGCCCGAGTGTGTAGGCCAGTTGCTGCAGGCCAGCTGCTGT
 GTGTCCAGCTGCTGCAGACCCAGTGTGCCAGTCTGCGTGTGCCAGCCACCTGCTGCCGCCCGAGT
 GCTGCCGCCCGAGTGTGCATCTCCAGCTGCTGCCAACCATCTTGTGGTGGTTCTAGCTGCTGTGGCTC
 CAGCTGCTGCCGCCCTTGTGCCGCCCTGTTGCCGCCCTGCTGCTGTCTGAGACCAGTCTGCCGTCAG
 GTCTGCTGCCAAACCACTTGCTACCGCCCACTTGTGTACCTCCACCTGCCCTGCCCCATGTGCTGTG
 CCACCCCTGTTGCT**GTA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_029613
Insert Size: 507 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_029613.1</u> , <u>NP_083889.1</u>
RefSeq Size:	978 bp
RefSeq ORF:	507 bp
Locus ID:	76444
Cytogenetics:	11 D