

Product datasheet for **MC211764**

1700123K08Rik (NM_029693) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 1700123K08Rik (NM_029693) Mouse Untagged Clone
Tag: Tag Free
Symbol: 1700123K08Rik
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC211764 representing NM_029693
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGTTCTCCTGTTTTAGGCTTCCAGAGGCTCCGGCCACAAAAAGCCAAGAGAGGTGCGCTTGTCTGCT
TTGGAGACGTTTGATTCGCCCTCTGACACACTTGAGGCATGCTTCCCACAGCGAGCCCAAGGTCTGTTG
TGAAAATGATCAGGAGCCAGACTCTATGCCAAGCACCCAGTCAACTACAAGAGTCGAGAATACGTT
CAACAGATGATCCACTACATACCAGCTGCTATTCAGAACCGTGACCACCTCTGTCTCGACATTTTTGTTG
CCATGTACCAAACGTATGCTACCACCTGGGAGGTGCTGGACCTGCTGATGAAAACGTATGCATCCTTCCA
GCCAGACTTCGTGGAAGACCAGCAAATAAGAGAGCCATATTCAGCTTTCTGTTCCGGCTGGTTTCAGAAA
TTCCCCCAGGATTTTTATGAGTCCCAGACCTGGCTGTTTTGAGTCAGTTCACGGAGTACGTGAGGCTCA
ATGTGCCTTCCCAGACATAGATACCCAAGCCAGGAGCCTCCCTCCATGTTGGAGGATCAGGAGGCCAT
CGCACTTAACTTGAGGAAGATCCCCTAGGGCAGGATGCCTCTGGGAGGCAGGAGACTCTGGCTCTGAGG
CCAGCTTCTATGGCAGAGCCACAGGGAGACGAAAAACCCGGAGAGAACACTGAGCTTGTGAAAAGAGCAG
TTCTGGATCCCTTGAACCAAAGCCACAATAGAGCTGCCTCCGTCTTACATCAAGCTCTGCCACATC
AGCTGATACATACTCACCTGTAGATGTGCAATAGATGAAGCTGCAGATGAGGCTGCAGATGTTTACCT
GCCAGTGTACAATCTACTTATGTTTTACCCGTGGACTGCAACTT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul
ACCN: NM_029693
Insert Size: 888 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_029693.2, NP_083969.1</u>
RefSeq Size:	1047 bp
RefSeq ORF:	888 bp
Locus ID:	76658
Cytogenetics:	5 G2