

## Product datasheet for MC211768

### Arpc2 (NM\_029711) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Arpc2 (NM\_029711) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Arpc2  
**Synonyms:** 34kDa; 2210023N03Rik; p34-Arc  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC211768 representing NM\_029711  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGATCCTGCTGGAGGTGAACAACCGCATCATCGAGGAGACGCTCGCGCTCAAGTTCGAGAACCGCGCCG  
CCGAAACAAACCAGAAGCTGTAGAAGTAACATTTGCAGATTTTGATGGCGTCTTTATCATATTTCAA  
TCCTAATGGAGACAAAACAAAAGTGATGGTCAGTATTTCTTTGAAATTCTACAAGGAACCTCAGGCACAT  
GGAGCTGATGAGTTGTTAAAGAGAGTGTATGGGAGCTTCTTGGTAAATCCAGAACCAGGATATAATGTTT  
CTTTGCTATATGACCTGGAAAATCTACCTGCATCCAAGGATTCTATTGTGCATCAGGCTGGCATGTTGAA  
GCGAAATTGTTTTGCCTCTGTGTTTGAGAAATACTTCCAATTCCAAGAAGAGGGCAAAGAAGGAGAGAAC  
AGAGCAGTTATCCATTATAGGGATGATGAGACCATGTATGTGGAATCTAAGAAAGACAGAGTCACAGTAG  
TCTTCAGCACGGTTTTTAAAGATGACGACGATGTGGTCATCGGAAAGGTGTTTCATGCAGGAGTTCAAAGA  
AGGACGAAGAGCCAGCCACACAGCCCCCAGGTCCTTTTCAGCCACAGGGAACCTCCTCTAGAGCTGAAA  
GACACAGATGCTGCTGTGGTGACAACATTGGCTACATCACCTTTGTGCTGTTCCCTCGCCACACCAATG  
CCACTGCTCGAGACAACCCATCAACCTGATCCACACGTTCCGGGACTATCTGCATACCCATTAAGTG  
CTCTAAGGCCTATATTCATACACGAATGCGAGCAAAAACATCCGACTTCTTAAGGTGCTCAACCTGCA  
CGCCAGATGCCGAGAAAAAAGAAATGAAAACAATCACGGGAAGACTTTTTTCATCCCGCTAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_029711  
**Insert Size:** 903 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_029711.1</a> , <a href="#">NP_083987.1</a>
<b>RefSeq Size:</b>	1597 bp
<b>RefSeq ORF:</b>	903 bp
<b>Locus ID:</b>	76709
<b>UniProt ID:</b>	<a href="#">Q9CVB6</a>
<b>Cytogenetics:</b>	1 C3
<b>Gene Summary:</b>	Actin-binding component of the Arp2/3 complex, a multiprotein complex that mediates actin polymerization upon stimulation by nucleation-promoting factor (NPF). The Arp2/3 complex mediates the formation of branched actin networks in the cytoplasm, providing the force for cell motility. Seems to contact the mother actin filament. In addition to its role in the cytoplasmic cytoskeleton, the Arp2/3 complex also promotes actin polymerization in the nucleus, thereby regulating gene transcription and repair of damaged DNA. The Arp2/3 complex promotes homologous recombination (HR) repair in response to DNA damage by promoting nuclear actin polymerization, leading to drive motility of double-strand breaks (DSBs).[UniProtKB/Swiss-Prot Function]