

Product datasheet for MC211797

Tsacc (NM_029801) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tsacc (NM_029801) Mouse Untagged Clone
Tag: Tag Free
Symbol: Tsacc
Synonyms: 1700021C14Rik; SIP
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC211797 representing NM_029801
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGAGCAGCACACTAGTAACCCTACTAAGGAAAAGCCAAAGAGGAAGATAGCGTTGCGTGTTTTTGTC
 CAGCAAAATCCTCCCCAGTTACATTGATCTTCCAGCAAATTTACGCCAGCTAATTTCTTGACCATCCA
 GACAAAGCTGTCTTCAGGAGCTGGTCAGAAGCCCAAGGGATGCCTCGGAGTCCTCGAGTGCATGTATGCC
 AACCTCCAGCTCCAGACCCAGCTTGCTCAACGACAGATAGCTATTTGGAAAACCTGCAGGCATCCCTCT
 CACAGCCGGGTTCTGGGAGGGAGAGCAAGAAGCTTCTCTCCCTGGCTTATGTTGCAAACCTGTTACTAAA
 TCACCTGCCCAATTCCATAAATGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: NM_029801
Insert Size: 375 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_029801.2](#), [NP_084077.1](#)

RefSeq Size: 481 bp

RefSeq ORF: 375 bp

Locus ID: 76927

UniProt ID: [Q9DA44](#)

Cytogenetics: 3 F1

Gene Summary: Co-chaperone that facilitates HSP-mediated activation of TSSK6.[UniProtKB/Swiss-Prot Function]