

Product datasheet for MC211798

Lypd5 (NM_029806) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lypd5 (NM_029806) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Lypd5
Synonyms:	2210003I03Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_029806, the custom clone sequence may differ by one or more nucleotides

```
ATGGGGATCCCCAGGACCTTCCTGCACTGCTTCCTCGGGAGCCTGCTCTGCCTGACAGGCTCCCAGGCC
TACAGTGCTACAGCTACGAGCACACCTACTTCGGGCCCTTCGACCTCAGCGCTGTGAACTCCCCAGTGT
CTCTTGTCTGAGGGGTGCTCTGAGGTGGTCTTGTGCTGGATACCGGTTACCGCAGCCTGGTGACAATG
GTCCGGAAAGGATGCTGGACAGGCCCCACCACCGGCCCATGCAGACCAACCAGGATGCTCTGCCTCCC
ACTACGGGTGGTCCGTGGTTCGCCACTGACTACTGCAACAGCAACCTCAAGACTCACGACGCCCTCCC
CAACCTGAGCCAAGCTCCTAACCACAGACGCTCAGTGGCACCAGTGTATGCCTGCTTGGGGACCCAC
CCTGAAGACTGTTCCCCTGAGAAGTCCCGACGGGTCCAGTGCCACCAGGACCAGAGCGCCTGTTTTCAGG
GCAACGGCCGGATGAACATTGGCAACTTCTCGGTGCCCGTGTACATCCGGACCTGTACCGGCCCTCCTG
CACCACCATGGGCACCACCAGCCCTGGACATCCATCGACCTGCAGGGTACTGTGCGAGGGCCACCTC
TGCAACAGGGCCTCGGTGACACAGATCCTGCCTGGCACCATGTCTCAGCCCTCCAGAGCTCCCCGAA
TCCTGACCCTGCTGACGGCTGCTCCACTACTGGCCATCGCTCTAGGGGCGTCCGTGCGGTTCCTGCATA
G
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_029806
Insert Size:	771 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC107186</u> , <u>AAI07187</u>
RefSeq Size:	935 bp
RefSeq ORF:	771 bp
Locus ID:	76942
UniProt ID:	<u>Q9D7Z7</u>
Cytogenetics:	7 A3