

Product datasheet for MC211808

Zfp934 (NM_175252) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zfp934 (NM_175252) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Zfp934
Synonyms:	6720457D02Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC211808 representing NM_175252 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGGAGACCTACAGGAACCTCAATGCTATAGGCTTTAATTGGGAAGCCAAAATATTGAAGAACATT
GTCAAAGTTCTAGAAGACACAGAAGGTGTGAAAGGAGTCAGAGTGCAGAGAAACCCTCTGAATACTCA
ACGTGATAAAGCCTTTGCTTTATCTGCTCACAGTCATGCCAAAGGCATGAAAGGATTCATACAGAAAAG
ATCCCCTCTGAAGATATCACTGTGTCGAAGACGTTTTACCTTACACTAGTCTCCAAGTACATAAAAGAA
CACAACTGGACAGAAACCCTATGAATGTAACAATGTGGTAAAGGATTTGCAAAGCCAGTCATCTTAA
AAGGCACGAAAGATATCACTGGAGAGCAACCCTTACATGTAATGAATGTGATAAAGTCTTTTCTTTT
AAACATCATCTCCAAACACATATAAGCATACGTAAGTGGAGAGAAACCCTTCAAATGTAATCTGTGTGATA
AATCCTTTTCAAACATAATCATCTTCAAAGTCATAGAAGAACACATACTGGAGAGAAACCCTACAAATG
TGATCAATGTGATAAAGCCTTTTCTGAGAAATGTAATCTCCAAACACATATAAGGATACATTCTGGAGAG
AAACCCTACAAATGTGATCAATGTGATAAAACCTTTTCTGAGAAATGTCATCTCCAAAACATATAAGGA
TACATACTGGAGAGAAACCCTACAAATGTAATCAATGTGATAAAGCCTTTTCTCAATACTACTCTTAA
AACTCATAGAAGAACACATACTGGAGAGAAACCCTACAAATGTAATCAATGCGATAAAGCCTTTTCTCAA
TACAGTAATCTTCAAACATAGAAGAACACATACTGGAGAGAAACCCTACAAATGTAATCAGTGTGATA
AAGCCTTTTCTATGAAATGTAATCTCCAAACCATGTAAGGATACATACTGGACAAAACCCCTTCAAATG
TAATCAATGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_175252
Insert Size:	990 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_175252.4</u> , <u>NP_780461.2</u>
RefSeq Size:	2103 bp
RefSeq ORF:	1326 bp
Locus ID:	77117
Cytogenetics:	13 B3