

Product datasheet for MC211832

Defb12 (NM_152802) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Defb12 (NM_152802) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Defb12
Synonyms:	9230103N16Rik; BD-12
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC211832 representing NM_152802

Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGAAGAATCTCCCTCAAACATGGCTCTGAGCAGAGAGGTGTTCTATTTGGTTTTGCTTTGTTCTTCA
 TTGTAGTTGAAGTCCATCAGGGTCCTGGGCAGGACTTGAGTATCCAGTCATTCCAGGAGGGGAGAT
 TGCTGTGTGCGAGACATGCAGGCTCGGCCGGGAAAAATGCAGGAGAACATGCATAGAGAGTGAGAAGATT
 GCGGGCTGGTGAAGCTGAACCTCTTCTGCTGTGAGAGAGGATCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_152802
Insert Size:	258 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_152802.3, NP_690015.1</u>
RefSeq Size:	997 bp
RefSeq ORF:	258 bp
Locus ID:	77674
UniProt ID:	<u>Q8K4N3</u>
Cytogenetics:	8 A1.3
Gene Summary:	Has antibacterial activity.[UniProtKB/Swiss-Prot Function]