

## Product datasheet for **MC211857**

### Samt1 (NM\_030036) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Samt1 (NM\_030036) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Samt1  
**Synonyms:** Gm15144; OTTMUSG00000019268  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC211857 representing NM\_030036  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGATCAACTCATCGTTGATACTCGGAGTGTCATCTACAGATTTATCAATGAGAAGGACTGCATCTACA  
AGCTGAGCGGCCCTCCTTAGCAGCTTAGCAGCTTTCATGCTTGAAGTACTTATTGCAAGCATCTTATCCTG  
GCGCCTATGGGAATTTGACAGCAATGTTGTGCAGTTTGTGTCCTTTGGCCTGTTTGAAGCATATTACCT  
CAGCAGTTAACATCTCTGGAACACTAACCAAGATGCTGGTATACACCCCTATTGATTCAACTTGAACA  
TTTCAACTGAATTTATGTATGCTCAGAACCTGGTAGTATGGGCCATTTTATGATGAAGCCTGTGTTCTGGT  
TTTCTGTGTAATGGCCATTAATAATCAGCTGCACGAAGAACCCTAGTGGAGATGCAGATATATTGCTAC  
AAGATCTCTGCCTTAATTTGTCTGTTAGCAGCATGTTTACATTTGTTTCTGTGATCTGGAACCACATGG  
TAGATTTTATGGCCACACCCTTTGACTTTCCATCTGATTTTCTGTTAAGAAAGAAGCTTTAACAAG  
CAAACATTTAACTGTGGTGTACCAGTAGGCCTCCTTAGCCACCATGTCACCTTTGGTGTGATCATG  
TTTCTCTGAGATAAGCGATTTGAAACTAAAGCGTCTGTGAAGGCCAATGATGCTTCCAAAATGGCC  
TTCTAGATGCGTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_030036  
**Insert Size:** 714 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_030036.1</a></u> , <u><a href="#">NP_084312.1</a></u>
<b>RefSeq Size:</b>	984 bp
<b>RefSeq ORF:</b>	714 bp
<b>Locus ID:</b>	78092
<b>Cytogenetics:</b>	X F3