

## Product datasheet for MC211874

### Lce6a (NM\_001166172) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Lce6a (NM_001166172) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Lce6a
Synonyms:	2210420J11Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC211874 representing NM_001166172 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTCACAGCAGAAGCCACAACCTTTCAGAGCTGCCAAATGCTCCCATATATTCACCTCCTAAAGGCCAA  
ACCTCTCTTTGAGTACCTGTTCAACTTCTTGTGGGACTTCTTGTTCAGCAGGCTGTTCTTCTACTCAA  
AAGGCTGAGGCTCCAGAACACTGTCAGCCATAAGGAAATTCATCACCCACAACCTCGCTGTCTTAGGGT  
AGTACCACCTACCACTGCAAAGAAGAAGAGTGC**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001166172
Insert Size:	246 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** NM\_001166172.1, NP\_001159644.1

**RefSeq Size:** 791 bp

**RefSeq ORF:** 246 bp

**Locus ID:** 78382

**Cytogenetics:** 3 F1