

Product datasheet for MC211882

1700057G04Rik (NM_001033184) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 1700057G04Rik (NM_001033184) Mouse Untagged Clone
Tag: Tag Free
Symbol: 1700057G04Rik
Synonyms: RP23-21J17.1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC211882 representing NM_001033184
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGCTATCAACACCACCTCATACGGAGTGCCACCTGGATTGGAATTTACTTCATATAGATCATATTA
 TGATTCAACAGCAGTTTGAGTTTGTGGAAGCCGTTTTAGGTTTTGAAACTGCTAACCGATATAAAATCAA
 TGACAAACTGGGACAGAAGGTTTACTATGCTGCAGAAGACTTTAATTTCTAACACTCAATTGCTGTGGG
 GCCATTAGACCTTTCACCATGAGGATCTTTGATAACTCAGGTCGTGAAGTCATAACTCTACGAAGACCAC
 TACGATGTGATTGTTGTTGCTGCCCTGCTGCCTCCAGCAGATAGAAATACAAGCTCCTCCTGGTGTGCC
 AGTGGGTATGTTATTCAAACCTGGCACCCATGTCGGCCAAAGTTTACAGTTCAAATGAGGAGAAACAG
 GATGTTCTAAAAATTATTGGTCCAATCTGCGTGTGCAACATTGGAGGAAGCATAGATTTTGAGATCAAAT
 CTCTTGATGAAGAATTCGTGGTTGGCAGGATTTCCAAGCACTGGTCTGGTATTCTGAAGGAGATACTGAC
 AGATGTAGACACTTTTGGGATCCAGTTCCTTTAGACCTCGATGTAAGATGAAGGCGGTGATGCTTGGT
 GCTTGTCTTCATAGACTTCATGTTTTTGAAGTGGTGGAGGCCAGAGACCAAACTTTTCTGG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001033184
Insert Size: 699 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001033184.3</u> , <u>NP_001028356.1</u>
RefSeq Size:	1052 bp
RefSeq ORF:	699 bp
Locus ID:	78459
Cytogenetics:	9 E3.3