

Product datasheet for MC211918

Klra23 (NM_024470) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Klra23 (NM_024470) Mouse Untagged Clone
Tag: Tag Free
Symbol: Klra23
Synonyms: Ly-49W1; Ly-49W2; Ly49W
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC211918 representing NM_024470
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGTGAGCAGGAGGTCACTTTCTCAGCTGTGAGATTCCATAAGTCTTCAGGGTTCAGAACCGGGTGA
GGCTTGAGGAGACAGGGAAGCCTCAAAAAGCTGGCCTCAGAGTGTGTTCCAGTTCCTGGCAGCTCATTGT
GATAGCTCTCGGAATCCTCATTTCCCTTCGACTGGTAATTGTCTCAGTGTGTTGACAAACATTTTCAG
AATAGTCAACAAAATCATGAACTGCAGGAACTCTAAACTGCCACGATAAGTGCAGCACCACCACGCAA
GTGACATCAACTGAAGGATGAACTGCTGAGTTCTACATCTATAGAGTGTAGGCCAGGCAATGATCTTCT
GGAATCCCTCCACAAGGAACAGAACAGATGGTACAGTGAACCAAGACTTTTTCAGATTCCTCACAGCAC
ACAGGCAGAGGTTTTGAAAAATATTGGTTCTGTTATGGTATAAAATGTTATTATTTTCGTCATGGACAGAA
AAACGTGGAGTGGATGTAACAGACCTGCCAGATTTCCAGCTTATCCCTTCTGAAGATAGACAATGAGGA
TGAAGTGAAGTTCCTTCAGAACCTGGCTCCTTCAGACATTTCTGGATTGGATTTTCATATGACAATAAG
AAAAAAGATTGGGTATGGATTGACAATGGCCATCTAAACTTGCCTTGAACACAACGAAATATAATATA
GAGATGGATTATGTATGTCGTTATCTAAAACAAGACTAGACAATGGGGACTGTGGTAAATCATACATCTG
TATTTGTGGTAAGAGACTGGATAAATTCCTCAT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_024470
Insert Size: 807 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_024470.1</u> , <u>NP_077790.1</u>
RefSeq Size:	942 bp
RefSeq ORF:	807 bp
Locus ID:	79410
Cytogenetics:	6 F3