

## Product datasheet for MC211948

### Vmn1r56 (NM\_030740) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Vmn1r56 (NM\_030740) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Vmn1r56  
**Synonyms:** V1rd3; V3R3  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC211948 representing NM\_030740  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGC**C

ATGGAAGTGTGGCTCTTCACATCCTTTTGCTTTGCCATGTTGTAGTTGGGACTGTGGCAACATCCTTC  
TGTTTGTCCATAATTTCTCTCCAATCTTGACTGATCCCGACTGAGGCCACATACAGGTCATTCTAATCAA  
CTTGGCTGTGGCAATGCATTATGCTCCTTCTCTTGACATATTCACATGACATGATTGATTTTGTCCA  
AAGAACCCTCAAAGTACCTCAAATGTAACTGCATACTTCTTTACATGGTGGCTCGAGGCACAAACA  
TGTGTTCTACCTGTGCTCCTCAGCACCTACCAGTTTGCCTCTTGTTCCTGGTAACTGCGCTAGGGTCAT  
GTTCAGAGAAATATCCCAAGGTCTGAGCTATTCTTGTTCAGTTGCTGGTTGTTCAAGTGTCTTAAAT  
AATGCTTATATCCCAATGAATGTGAGTGGTCCACAGAAATCACACAATGACACTGATTCTAAAGGCAATT  
GTATATGCTCCGTCTCTGGGTCAGTGTAGACATGAATTTCTTTTCGTTTTCCCATGACATCATGTTCT  
TAGCATCATGGCTGGACCAGTGTCTCCATGGTATTCAATTAACAGACACCACCAGAGAATGCACTAT  
ATACATAATCATAACCATAAAAAACAGGCCATGCTGAAACCAGAGCAGCCACACCATCTGATGCTGG  
TGGTCACATTTGTGAGCTTATATATTCTAGATTGTATTTGATTTTATTTACATTTCTTTGTGGAAATC  
TCGACTCTGGTTGAGGCATGTCATAAAGTATGACTCTAAGCTTCCCCACCATTCTCCCTTACTGTTG  
ATTTTTAGGGATCCTAGGGTCTATTGTTCTCTGCTCTTCAAGTGTGGGCTGGAAATCCACATGACTAGGG  
GAGACATTACAGAGGAACAG**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_030740  
**Insert Size:** 933 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_030740.1</a> , <a href="#">NP_109665.1</a>
<b>RefSeq Size:</b>	1836 bp
<b>RefSeq ORF:</b>	933 bp
<b>Locus ID:</b>	81015
<b>Cytogenetics:</b>	7 A1