

## Product datasheet for **MC211992**

### **Uchl4 (NM\_033607) Mouse Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Uchl4 (NM\_033607) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Uchl4  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC211992 representing NM\_033607  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGAGGGTCAACGCTGGCTGCCACTGGAGGCCAACCCCTGAGGTCACCAACCAGTTTCTCAAGCAATTAG  
GCCTGCATCCTAACTGGCAGTTTGTGATGTACGGGATGGAGTCTGAACCTTCTAGCATCATACCAAG  
ACCAGTGTGCGCAGTGTACTCCTCTCCCTATCACAGAAAAGTATGAAGTCTTCAGAACAGAAGAGGAA  
GAAAAGATAAAATCTCAAGGACAAGATGTGACATCATCAGTATATTTTATGAAACAAACCATCAGCAATG  
CCTGTGGAACGATTGGAACGATTGGACTAATCCATGCCATTGCAAACAACAAAGACAAGGTGCACCTTTGA  
ATCAGGGTCAACATTGAAAAAGTTCTGGAGGAGTCTGTATCAATGAGCCCTGAAGAGAGAGCCAAATAC  
CTGGAGAACTATGACGCTATTCGAGTTACTCATGAAACCAGTGCACATGAAGGTCAGACTGAGGCACCAA  
GTATAGATGAAAAAGTAGATCTTCATTTTATTGCGTTAGTACATGTAGATGGGCATCTCTATGAATTAGA  
TGGATGGAACCATTTCCAATTAACCATGGGAAAACCTAGCGATGAGACGTTGTTAGAGGATGTCATAAAA  
GTTTGCAAGAAGTTCATGGAACGTGACCCTGATGAGTTAAGATTTAATGCAATTGCTCTCTCGGCAGCAT  
AG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-MluI  
**ACCN:** NM\_033607  
**Insert Size:** 702 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_033607.1</a> , <a href="#">NP_291085.1</a>
<b>RefSeq Size:</b>	1162 bp
<b>RefSeq ORF:</b>	702 bp
<b>Locus ID:</b>	93841
<b>UniProt ID:</b>	<a href="#">P58321</a>
<b>Cytogenetics:</b>	9 34.58 cM
<b>Gene Summary:</b>	Ubiquitin-protein hydrolase is involved both in the processing of ubiquitin precursors and of ubiquitinated proteins. This enzyme is a thiol protease that recognizes and hydrolyzes a peptide bond at the C-terminal glycine of ubiquitin.[UniProtKB/Swiss-Prot Function]