

Product datasheet for **MC212001**

Tm2d1 (NM_053157) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tm2d1 (NM_053157) Mouse Untagged Clone
Tag: Tag Free
Symbol: Tm2d1
Synonyms: 2310026L18Rik; AA990549; Bbp
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC212001 representing NM_053157
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGCCGCTGGCCCGGGTCGGGCTTCCCAGCGGGGGCTCCGGGCTTCTCCGCACCCTGT
GGCTCGTGACGGTCGCCCGGGACACTGTGGGCTGCTGCCTCTGGCGCTGTCGGGGCGAGGAGACACC
CAAGTGTGAGGACCTCAGGGTGGACAATATTTGTAAGAACCAAAAATAAATGATGCTACGCAAGAA
CCAGTTAATTGTACAACTACACAGCTCATGTTCAATGTTTTCCAGCACCAAAAATAACTTGAAGGATT
TGAGTGGTAATGAAACACATTTTACTGGAAGTGAAGTCGGTTTTCTCAAGCCCATATCTTGCCGAAATGT
GAATGGCTATTTCGTACAAAGTGGCAGTTGCATTATCTCTTTTTGGGATGGCTGGGAGCAGATCGATT
TACCTCGGATATCCTGCCTTAGGCTTGTAAAAATTTGCACCGTAGGATTTTGGGAATTTGGGAGCCTAA
TTGATTTCACTTATTTCATGAGATTGTTGGACCTTCAGATGGAAGTAGTTACATTATAGACTATTA
TGAACCAGGCTTACAAGACTCAGCATTACTAATGAAACATTTAGAAAAACCCAGCTGTACCCAT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: NM_053157
Insert Size: 627 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_053157.2</u> , <u>NP_444387.1</u>
RefSeq Size:	1080 bp
RefSeq ORF:	627 bp
Locus ID:	94043
UniProt ID:	<u>Q99MB3</u>
Cytogenetics:	4 C6
Gene Summary:	May participate in amyloid-beta-induced apoptosis via its interaction with beta-APP42. [UniProtKB/Swiss-Prot Function]