

Product datasheet for MC212003

Mrp11 (NM_053158) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	Mrp11
Synonyms:	2410002L03Rik; 5830418D04Rik; AI462572; Llmt
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC212003 representing NM_053158 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCGGGCGCCGTAAGGTGCCTCCGTAGAGTCTTGATACATCATCAAAGGCATTGTCTTTGCAAGATGG
CTTCCCAGGCATCTCTTTACCTTGTCTGTGAACAGTCTCTTGCACAACAGACATTTTGCTGCTGCTGC
TGCTGCTGTACAAAGCCTGCAAGAAAAATTAAGGAGTCCCAAGGAAAGACATCAGATGAAAAACCA
GTAGATGATATAGAGAAAATAAGTCATATACGTACATGGAAGTGATCCTGAAGATGATGTCTACTTAA
AACGCTTATACCCACGGCGGATCTATGAAGTAGAGAAAGCTATTCACCTTACTTAAGAAGTTTCAAGTTCT
GGACTTTACCAATCCAAGCAAGGTGTTTATCTTGACTTAACATTAGATATGGCATTGGGGAAAAAGAAA
ACAGTGGAGCCATTTGCCAGTGTTATTGCTTTACCACACCTCTTCTCCTCAGAAGTCAACAAGTTGCTG
TGTTTACAGCGAACGCATCAGAAATTAAGTAGCAGAGAAATGGAGCTGCGTTTGCAGGAGGAACAGA
TCTGGTTAAGAAGATCATGGATGACGAGGTTGTTGTGGACTTTTATGTAGCTGTCCAGAAATAATGGGT
GAACTTAACCCCTTAAGGAAGAACTGAAAAAAGATTTCCAAAGGCTACTAGAAATCCATTGGTCGTG
ACATCCCCAAATGCTTGAATTATTTAAACGGCACATGAAATTATGGTAGATGAAGAAAGACAGAACTT
CCTCAGCACCAAAATAGCAACATTGGATATGCCAAGTGACCAGATAGCTGCCAATCTTCAAGCAGTTATA
AATGAGGTCTGTAAAGCACAGACCACTGAATCTGGGTCCGTTTGTGGTCCGTGCTTCTCTCGTAGTTCAA
CCAGTGAAGGTTTACTACTGAAGACTGACTCATTGCTACCCAAAGAAGCTAAACTACAGAAGCTGAAAC
GGAAGAACTCAACGGCAGAAGCTGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



ACCN:	NM_053158
Insert Size:	1011 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_053158.3</u> , <u>NP_444388.2</u>
RefSeq Size:	3449 bp
RefSeq ORF:	1011 bp
Locus ID:	94061
UniProt ID:	<u>Q99N96</u>
Cytogenetics:	5 E3