

Product datasheet for **MC212026**

Fam8a1 (NM_001033192) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fam8a1 (NM_001033192) Mouse Untagged Clone
Tag: Tag Free
Symbol: Fam8a1
Synonyms: C78339
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC212026 representing NM_001033192
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCAGAGGGGCCGAGGAAGCCCGGAGCCGCCCTCCGGAGCAGGATGGCAGTGGCGGGGACCACGAGC
CGGTCCCTTCGAGCCAGGCTCTCCACCACCACCGTCGCGGCCACCGCCGCCCCCCCGGTCCCCGGC
CCTCGCGAGCGTGCCTCGGCAGTGTGGAGCCCGAGCCGAGTCCGAGCTCCGCAAGCGCCGGGAGACG
GGCGCGGAGCCTCCGGAGGCGGGAGACTGCGAGGCGTCCGACGGCCCGGGGCGGCTGAGCGCGCGGAGT
ACTCGCGCAGGTGCACGAATGGCTGTGGCAGTCTACTGCGGCTACCTGACCTGGCACAGCGGCCCTGGC
CGCCCTGCCCGCCTACTGCGGCCCGCCACCACCCGCGCCACCGCCCGCAGCCCGCCGTCCCCTCGCCG
TCGTCCCCGCGGACGCGCTCGTCCCCGCGCCACCGCGCCGCGCCAGCTGGCGTACTACAACCCCTTCT
ACTTCTGAGCGCTGCGGGGCCGGGGCGGGGCGGCCACGGGCGGCGCCACCCGACGGCTGTGCGGGG
CCTCACGGCCCGGGCCCCGACGTGCAGCCATCGGCCGGGCGGCCCGGGGACCAGGGTGGGTTCTGCA
GCCCCGCGCGGACTGCCAGCGACACCGGGCGGCGAGGAGTATGATGATGATGATGATGATGATGATGATG
ACAGATTTATGGCAGAGATGGTGGATTTCTTTATTCTCTTTTATAAAAGCAACCATTTGCTTTGAGCAT
CATGCATCTCAGTGAATAAAGGACATCTCTAAGTTTGCCATGCACTATATAATAGAAGAGATTGATGAA
GACACATCAATGGAAGACCTGCAGAAAATGATGATTGTGGCGCTTATATACAGGTTGTTAGTTTGTCTT
ATGAGATAATTTGCATTTGGGGAGCAGGGGAGCTACCCAGGAAAGTTCTGCTAGGGCTTCGAGTTGT
GACATGTGATACATCAGTACTCATTGCTCCAAGTCGGGTTTTAGTAATTCCTTCTCAAATGTTAGCATC
ACTACGTCACCATACGAGCTTTGATCAAGAATTTTTCTATTGCTTCTTTTTTCCCTGCTTTCATCACAC
TGCTGTTTTTTCAGCATAATCGAACAGCCTATGACATTGTAGCAGGAACCATTGTGGTAAAGAGAAATGG
GGTCAGATGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001033192
Insert Size:	1200 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001033192.2</u> , <u>NP_001028364.1</u>
RefSeq Size:	1381 bp
RefSeq ORF:	1200 bp
Locus ID:	97863
Cytogenetics:	13 A5