

Product datasheet for **MC212124**

Ypel1 (NM_023249) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ypel1 (NM_023249) Mouse Untagged Clone
Tag: Tag Free
Symbol: Ypel1
Synonyms: 0610009L05Rik; 1700016N17Rik; 1700019O22Rik; 4921520K19Rik; 4930511F14Rik; AV118478
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC212124 representing NM_023249
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGTGAAATGACCAAGTCCAAACTTTCCAAGCTTATCTGCCACACTGTCACCGAACTTACAGCTGCA
TCCACTGCAGAGCCCATCTGGCCAATCAGGATGAGCTTATCTCCAAGCTTTTCAGGGAAGCCAGGGGCG
AGCCTACCTCTCAACTCTGTGGTGAATGTGGGCTGTGGCCCGCTGAGGAGAGAGTCCTTCTAACC
CTGCATGCAGTGGCCGACATCTACTGTGAGAACTGCAAGACCACACTCGGGTGGAAATACGAACATGCCT
TTGAGAGCAGTCAGAAATACAAGGAAGGAAAATTTATTATTGAACCTGCCACATGATCAAGACAACGG
CTGGTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: NM_023249
Insert Size: 357 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_023249.6](#), [NP_075738.1](#)

RefSeq Size: 2879 bp

RefSeq ORF: 357 bp

Locus ID: 106369

UniProt ID: [Q9ESC7](#)

Cytogenetics: 16 A3

Gene Summary: May play a role in epithelioid conversion of fibroblasts.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (2) lacks an alternate in-frame exon in the 5' coding region and uses a downstream start codon compared to variant 1. It encodes isoform 2 which has a shorter N-terminus compared to isoform 1.