

## Product datasheet for MC212162

### Fxyd4 (NM\_033648) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Fxyd4 (NM\_033648) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Fxyd4  
**Synonyms:** 0610008I02Rik; A1267073; Chi; Chif  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC212162 representing NM\_033648  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGAGGAAATCACCTGTGCCTTCTCCTGCTGCTAGCAGGTCTGCCGGCCTTGAAGCCAGTGACCCAG  
TTGATAAAGACAGTCCCTTCTACTATGACTGGGAGAGCCTGCAGCTGGGAGGATTGATTTTTGGAGGGCT  
CCTGTGCATCGCTGGAATTGCCATGGCCCTGAGTGGCAAGTGCAAATGTAGGCGTACCCATAAGCCAGT  
TCCTTACCTGAAAAGCCACTCCACTCATCATTCCAGGCTCTGCCAATACCTG**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-MluI  
**ACCN:** NM\_033648  
**Insert Size:** 267 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_033648.1](#), [NP\\_387468.1](#)

**RefSeq Size:** 554 bp

**RefSeq ORF:** 267 bp

**Locus ID:** 108017

**UniProt ID:** [Q9D2W0](#)

**Cytogenetics:** 6 F1

**Gene Summary:** This gene encodes a member of a family of small membrane proteins that share a 35-amino acid signature sequence domain, beginning with the sequence PFXYD and containing 7 invariant and 6 highly conserved amino acids. The approved human gene nomenclature for the family is FXYD-domain containing ion transport regulator. Mouse FXYD5 has been termed RIC (Related to Ion Channel). FXYD2, also known as the gamma subunit of the Na,K-ATPase, regulates the properties of that enzyme. FXYD1 (phospholemman), FXYD2 (gamma), FXYD3 (MAT-8), FXYD4 (CHIF), and FXYD5 (RIC) have been shown to induce channel activity in experimental expression systems. Transmembrane topology has been established for two family members (FXYD1 and FXYD2), with the N-terminus extracellular and the C-terminus on the cytoplasmic side of the membrane. [provided by RefSeq, Jul 2008]  
Transcript Variant: This variant (1) represents the longer transcript. Variants 1 and 2 encode the same protein.