

Product datasheet for MC212171

B3gnt5 (NM_001159407) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	B3gnt5 (NM_001159407) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	B3gnt5
Synonyms:	beta3GnT5
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>MC212171 representing NM_001159407 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAGACTGTTTGTAGCAGAAGAGTCAAAGATGGAAAATTTTCACTTTTTGTCACCTGTTTTATAT
TAAGCTTCATGGTTTTTGGAGCCCAATCAATAATTACATCATGAGCCATATGAAGTCTACTCCTACAG
ATACCTCGTAAATAGCTATGGCTTTGTAACAATTCCTGTCTCTCAAGCACAGCTCTGTGCAGCCTCAC
TACCCATACTTGATCAACCACAGAGAGAAGTGTCAGGCTCAAGATGTCCTCCTTACTGTTATAAAGA
CTGCCCTGAAAACATGCGCCGACGTTCTGCAATCAGAAAAACGTGGGGCAATGAGAACTATGTTACGTC
TCAACTCAATGCCAACATCAAAATCTGTTTGCAATAGGAACTCCTGGTCCACTGAAGGGAAAAGAAGT
CAAAAAAGACTAATCGGGGAAGATCAAGTGTACAAGGATATAATTCAGCAAGATTTCAATTGATTCTTTCC
ACAATCTTACTTCTAAATTCCTTCTTCAAGTTCAGCTGGGCAAATACCTTTTGTCCACATGCCAAATTCCT
GATGACTGCTGATGATGATATTTATCCACATGCCAAATCTCATTGAATATCTCAAGGGCTAGAGCAG
ATTGGAGTTCGAGACTTTTGGATTGGTCACGTTTCATCGAGGTGGCCCTCCTGTTAGGGATAAAAGCAGCA
AATACTATGTTCCCTATGAAATGTACAAGTGGCCAGCCTACCCTGACTATACAGCTGGTGCCTATGT
TGTCTCCAGAGATGTAGCTGCCAAATCTATGAGGCATCGCAGACTGAATTCAGTATGTACATAGAT
GATGTATTCATGGGCCTCTGTGCCAATAAAGTGGGAATCTTGCCACAGGACCATGTATTTTCTCTGGGG
AAGGGAAAATTCCTTATCACCCCTGCATCTATGAAAAGATGATGACATCTCACGGACACTTACAAGATCT
GCAGGACCTCTGGATAGAGGCTACACATCTAAAGTAAAGAACATTTCAAAGGGTTTTTGGTCAAATA
TACTGCAGGTTAATTAAGATAGTTCTTCTCTGCAGACTGACTTACAGGAATTCATACCCTGTTGGGCTG
CATTTGCT**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul



[View online »](#)

ACCN:	NM_001159407
Insert Size:	1131 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001159407.1</u> , <u>NP_001152879.1</u>
RefSeq Size:	4872 bp
RefSeq ORF:	1131 bp
Locus ID:	108105
UniProt ID:	<u>Q8BGY6</u>
Cytogenetics:	16 A3
Gene Summary:	<p>Beta-1,3-N-acetylglucosaminyltransferase that plays a key role in the synthesis of lacto- or neolacto-series carbohydrate chains on glycolipids, notably by participating in biosynthesis of HNK-1 and Lewis X carbohydrate structures. Has strong activity toward lactosylceramide (LacCer) and neolactotetraosylceramide (nLc(4)Cer; paragloboside), resulting in the synthesis of Lc(3)Cer and neolactopentaosylceramide (nLc(5)Cer), respectively. Plays a central role in regulating neolacto-series glycolipid synthesis during embryonic development.</p> <p>[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longest transcript. All three variants encode the same protein.</p>