

## Product datasheet for MC212177

### Zdhhc15 (NM\_175358) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Zdhhc15 (NM\_175358) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Zdhhc15  
**Synonyms:** 6030457O13Rik  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC212177 representing NM\_175358  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGCGCGGGGCTGGAAGATGGCTCTGTCTGGGGGTTGCGGTGCTGCCCGGGTACTGTCCTGGGTGC  
CAGTGCTCGTTATTGCCTCGTCGTTCTCTGGTCTACTATGCCTACGTCTTTGAACTCTGCCTGGTGAC  
TGTTTTGAGTCCAGCGGAAAAAGTTATCTACCTCATACTCTATCATGCCATTTTTGTGTTCTTTGCATGG  
ACCTACTGGAAGTCTATCTTCACTCCACAGCAGCCAAACCAGAAGTTCCACTTGTCCATACAGACA  
AGGAGCGCTATAAAAAATGAAGAGAGACCTGAGGTCCAGAAGCAGATGCTTGTGGACATGGCCAAGAAGTT  
ACCAGTGTACACAAGAACCAGGAGTGGAGCTGTACGGTCTGTGACCGGTGCCATCTAATCAAGCCAGAC  
CGCTGCCATCACTGTTCTGTCTGTGCCATGTGTGTTTTAAAAATGGACCATCACTGCCCATGGGTTAATA  
ACTGCATTGGATTTTCCAACACAAATCTTCCCTCAGTCTTAGCTTATCCGTTCTCTACTGCCTGTA  
CATTGCTACAACAGTCTTCAGCTACTTCATCAAATACTGGAGAGGGGAATTGCCAGTGTTCGCTCTAAG  
TTCCATGTCCTCTTTCTCCTCTTTGTGGCTGCATGTTTTTTGTCAGCCTGTGATTCTCTTTGGTTACC  
ATTGTTGGCTTGTGAGCAGAAACAAACTACCTTAGAGCCTTCTGCACACCAAGTATTCACAAGTGGACC  
AGAGAAAAATGGGTTCAATCTTGGCTTCATCAAGAATATTCAGCAGGTGTTTGGCGATAATAAGAAGTTC  
TGGTTAATACCCATTGGGTCCAGCCCTGGTGATGGACACTCCTTCCCTATGAGGTCTATGAATGAGTCAC  
AGAACCCGCTGCTGGCAAATGAGGAACCTGGGAAGACAATGAGGATGATAGTCGAGATTATCCAGAAGG  
CTCATCATCACTTGTGTGGAGTCGAAACGT**AG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_175358



[View online »](#)

<b>Insert Size:</b>	1014 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_175358.4</a> , <a href="#">NP_780567.1</a>
<b>RefSeq Size:</b>	5517 bp
<b>RefSeq ORF:</b>	1014 bp
<b>Locus ID:</b>	108672
<b>UniProt ID:</b>	<a href="#">Q8BGJ0</a>
<b>Cytogenetics:</b>	X D
<b>Gene Summary:</b>	Catalyzes palmitoylation of Cys residues on target proteins. Catalyzes palmitoylation of GAP43 and DLG4/PSD95.[UniProtKB/Swiss-Prot Function]