

## Product datasheet for **MC212232**

### Cebpe (NM\_207131) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cebpe (NM\_207131) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Cebpe  
**Synonyms:** C/EBPe; CRP1; Gm294  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC212232 representing NM\_207131  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGTCCCACGGGACCTACTATGAGTGCGAGCCTCGGGTGGCCAGCAGCCACTTGAGTTCTCAGGGGTC  
GAGCTGGGCCCGGGGAGCTGGGGACATGTGTGAGCATGAGGCCATCCATCGACCTCTCCGCTACATCGA  
GTCTGGGAAGAAGCAGCTACTTTCTGACCTCTTTGCCATGAAGCCAACGCCTGAAGCCCGAAGCCTAAG  
GGCCAGGAGCCCTCGTTTCCCTCACTACCTGCCGGTGAACCTCGGCCATTTCGCTATCCCTCACACA  
CATTTGGCCAGATAGGAAGGCTTTGGGCCTGGCATCTACAGCAACCCAGGGAGCTACGACCCAGGGC  
TGTGGCGGTGAAGGAGGAGCCTCGTGGCCAGAGGGCAACCGAGGCACCAAGTTCGAGGCAGCTACAATCCC  
CTGCAGTACCAAGTGGCACACTGCGGGCAGACAGCCGTGCACCTCCCCCGACCTGGCAGCACCTGGCC  
AGCCCTTGGTGTCTCAAGGCCCTGTGGCAGCTGCTGCCCCCTCCCTGCAGTCCCCTTCTCAAGGCACC  
CTCCCCGCTGGCCCCTCACACAAGGGCAAGAAGGCAGTGAACAAAGATAGCCTGGAGTACCGACTGCGA  
CGTGAACGTAACAACATCGCGGTGCGCAAGAGCCGGGACAAGGCCAAGAGGCGCATTATGGAGACTCAGC  
AGAAGGTCTGGAGTACATGGCTGAAAACGAGCGTCTGCGCAACCGCGTAGACCAGCTCACTCAGGAGCT  
AGACACTCTGCGAACCTCTCCGCCAGATCCCTGAGGCTGCTAGCCTCATCAAGGGTGTGGGGGCTGC  
AGCTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_207131  
**Insert Size:** 846 bp



[View online »](#)

<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>Components:</b>	<p>The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).</p>
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<p><a href="#">NM_207131.1</a>, <a href="#">NP_997014.1</a></p>
<b>RefSeq Size:</b>	<p>1238 bp</p>
<b>RefSeq ORF:</b>	<p>846 bp</p>
<b>Locus ID:</b>	<p>110794</p>
<b>UniProt ID:</b>	<p><a href="#">Q6PZD9</a></p>
<b>Cytogenetics:</b>	<p>14 27.9 cM</p>
<b>Gene Summary:</b>	<p>Transcriptional activator. C/EBP are DNA-binding proteins that recognize two different motifs: the CCAAT homology common to many promoters and the enhanced core homology common to many enhancers. Required for the promyelocyte-myelocyte transition in myeloid differentiation.[UniProtKB/Swiss-Prot Function]</p>