

Product datasheet for MC212254

Vmn1r44 (NM_053227) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Vmn1r44 (NM_053227) Mouse Untagged Clone
Tag: Tag Free
Symbol: Vmn1r44
Synonyms: V1RA10; V1rb4; V1rb11
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC212254 representing NM_053227
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAATAAAGCCAACCTTACTCCACATTGACACAAACATAAAAAACACCTTGCTTGCTGAAGTGAGTGTTG
GGATCTCAGCTAACAGTATCCTTTTATAGCCTACCTCTGCATGCTCCTGGTGAGAACAGGCATAAGCC
CATTGATCTCTACATTGCTTTCTTGTCCTAACTCAACTAATGCTGCTTATAACTATGGGACTCATAGCT
GTGGACATGTTTATGCCTTGGGGGAGATGGGATTCTACCACATGCCAATCCCTTATCTATTTGCACAGGT
TTCTGAGGGGCTCACCTTTGTGCTACTTGTCTGCTGAATGTCCTTTGGACCATCACCTTAGTTCTAG
AAACTCTTGTTTAGCAAAGTTTAAACATAAATATCCCCATCACATCTCAGGTGCCTTTCTTTTCTCTGT
GTTCTCTATATGTCTTTCAGCAGTCACCTTTTAGTATCAATGACTGTGACCCCAATTTGACCTCAGAGA
ATTTTATGTATGTTACTCAGTCCTGCTCACTTCTACCCATGAGTTACTCCAGAACAGCATGTTTCCAC
ACCGGTAGCCATCAGGGAACCTTTCTTATCAGTCTCATGGCCCTGTCCAGTGGGTACATGGTGGCTCTC
CTATGGAGGCACAAGAAGCAGGCCAGCATCTTCGCAGCACCAGCCTTTCTTCAAAGCATCCCCAGAGC
AAAGGGCCACCAGGACCATCTGCTGCTCATGAGTTTCTTTGTGGTTTCTACATTTGGACACTGTTAT
CTTCCACTCAAGGATGAAGTTCAAGGATGGCTCAATATTGTACTGTTTCCAAATTATTGTGTCCCATAGC
TATGTCACAGTCAGCCCTTTTGTGTTTATTGCACTGAAAAGCATATAATTAAGTTTGTGAGGTCAATGT
GTGGCAGAATAGCAAATATTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_053227



[View online »](#)

Insert Size:	933 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_053227.2 , NP_444457.2
RefSeq Size:	933 bp
RefSeq ORF:	933 bp
Locus ID:	113854
UniProt ID:	Q9EQ47
Cytogenetics:	6 D1
Gene Summary:	Putative pheromone receptor implicated in the regulation of social and reproductive behavior.[UniProtKB/Swiss-Prot Function]