

Product datasheet for MC212262

Vmn1r15 (NM_053236) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Vmn1r15 (NM_053236) Mouse Untagged Clone
Tag: Tag Free
Symbol: Vmn1r15
Synonyms: V1rc6
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC212262 representing NM_053236
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCCTCGGTAGAGAATGTCCTTTATTTCCAAGCTGGGTTTGGAGTCTACCCAATGTGTTTCTTCTTC
 TTTTCTATATTTTCATAATCCTAGGTCACAGACCTAAACCCACAGACCTGATCTCCTGTCAACTGACTTT
 GGTTACATAATCATGTTCTCCTCACTGGAGGAAATATTTGGCTCGCAGACATATTTGGGTTTCTAAACATT
 GAGAATGACATAAAATGTAAGCAACTTTTATACAAGTAGGGTAACGAGAGGCCCTCTATCTGCATTA
 CCTGCCTCCTGAGTGTTCAGGCTGTACAATCAGTCCAGTACCTCTATGCTGGCAAAATTTAAACA
 GAAACTTAAAAAGTACATGAACAATGCATTGTTCTATATTTGGTTTTCAATTTGTCTGTTGTAGTAAC
 ATGCTTTTCTTTGTTGGTGGTTTTACCAATGTCAGTGAGACCAAGCAGGTGAAAGTCACTAAATCCTGCT
 CACTCTTCCCATGAACATACATCAGGGAATTAATTGTTATTCTGTCAACCTCCAGAGATGTGTTTCT
 TGTAGCAGTCATGCTGACCACAAGTGCATATATGGTGAATATTTGCACAGGCATCAGAGACAATACAAA
 CATCTTCATAGCATCAAAGCCTGAGAACATCCCCTGAGAAAAAGGCCACCCAGACCATTGTGCTAG
 TAGTTTTCTTTGTAGTCATGTAAGTGGGATTTTCATCTCATTACCTCACTCTGTTATGATGTA
 TGACCCAGTCGTCCTGACAGTTCAGAAGTTTCTGATATATGCTTATCCCAAAATTAATCCTTTGGTACAA
 ATTAGTCTGATAAGAGAATAATAAATGTGCTGAAAAATTTGCAGTCAAAGTGCCACT**AG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_053236
Insert Size: 900 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_053236.2</u> , <u>NP_444466.2</u>
RefSeq Size:	900 bp
RefSeq ORF:	900 bp
Locus ID:	113863
Cytogenetics:	6 B3