

Product datasheet for MC212376

Vmn1r232 (NM_134193) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Vmn1r232 (NM_134193) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Vmn1r232
Synonyms:	V1re4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC212376 representing NM_134193 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGCTGTGTGACTTGAAATAGGTAGTACTGTGATTCTTGCACATACATGTCCAACTTAGAAACTTAC
TTCAATACATAGTTCACAAATGTGATTATTACCATAATTTAAAGGAAAAACGTTGAATGACCTAATGGA
GTCCAGGAACATGGCAATAGGAGTAGGACTCTCACTTCAGAGTATACTTGGTATTCTGGGAAATTTCTCT
CTTCTTTTCTACTACCTAATCCTATACTACAAAGAACACACATTAAGATCAGAGACATGATTCTAGTAC
ATGTGCTCACATCTAACTCCTTGATCATTATCTCAAAGGTTTCCCCAGATAATGGGAGCTTTTGGGTG
GAACAGTTATTCAATGATGTTGGGTGCAAACTTATTTTATATGTTCAAAGACTTAGCAGGAACATGTCC
ATTATCACCACCTTGCTCTTGTAGTGTCTCCAGGCAATCACCATCAGCCCCAGGAATTCCTATTGGAAGA
AATTTAAAGTCAAATCAACAAAGTTTATGGGCTTGCCCATTTCTCTCTTCTGGATCTTGTATTGTAGT
AAATATGCTTTTCCCTGTGTATACATCTACCAAGACGAATAGCAAAAACAGACAAAAAGAGAGATTCT
GAATTCTGCCATTCTCTAGGCCGTGACAAAATAGTTGATATACTATACGGCACTTTGTGTGTTTCTCTG
AAGTTGTATTTTCTCTGCTCATTGTATGTTCTAGCACCTTCATGATTGCAATACTTCATGTACACAAAA
GAGGGTTCAACACATCCTCTGTACTCATACTTATCCAAGAAGATCCCCTGAAAACAGAGCCACACAGACC
ATCTTGATTCTTGTGTTGACCTTTCTAGCATTTTATACACTCTCTCAGTTTGTCAAGGCTATATTGCTC
TTTCACATGATCCAGTTGGTGGGTAATGAATACCACAGCCATCTTTCTATGGGTTTCTCTACTTTAGG
GCCATTTGTGATGAGTCGTGATTCTACTGTTTCCAGATTCTGCTTTACCTGTATAAGAAATATAGAGAAA
ACGTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_134193
Insert Size:	1056 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_134193.2, NP_598954.2</u>
RefSeq Size:	1159 bp
RefSeq ORF:	1056 bp
Locus ID:	171227
Cytogenetics:	17 A3.2