

Product datasheet for MC212396

Vmn1r199 (NM_134213) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids

Tag: Tag Free

Symbol: Vmn1r199

Synonyms: Vlrh4

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC212396 representing NM_134213
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGCATCCACCACTGCCTCTTGAACACCAGAAAATGCTCTCATGCAAAATGTGCTGTTCCAGATCCATGTC
CAGCATTACTGATTCTTTCTTTAAAATGATTTGCATCCTGCACACAAGGAAGGAAGAAATGATTTTGAA
GTTAATTAAGAAAATACTTTTCTCTTCATGACTGTGCTTGGCATTATGGGAACATGTCTGTTTCTGTG
AATTACATGTTCACTGTTGTTGGGAAGCCCTGAGAAGAAACCCGTACACCTTATTCTTATCCACCTGGCTT
TTACAAACACTGTCACTCTTCTTGCAAAAGGACTGCCAAAGACTATTGTAGCTTTTGGTTTGAGAAATTT
CCTAGATGACTCAGGCTGTAAGATCATTGTTTACTTGAAAGGGTGGCCCATGGCCTCTCTCTGCGCC
AGCAGTCTCCTCACTGTGGTCCAGGCCATCATCATCAGCCCCAGAGCATCTGGATGGAGGATGCTCAGGC
CAAAGTCTGCATGGCACATCCTTCCATTCTTCTTATTCTTTGGTTATTCACTGCATTAATAAATATGAA
CCTAATAAATTCATCACAAGTACAAGACTGAATATATCACAATTAAGACTGATAGCAATAATTGTTAT
TTTATGTTACCAAGTCAGAAAATAAAGGGATTGTTCTCCCTCTCATGGTCTGAGAGATGCAGTCTTTC
AGGGTGCTATGGGAGGGGCCAGTGGCTACATGGTATTTCTTCTCTATAAGCACCACCAGCATGTCCTCTA
CCTTCAGAAGTCCAAGCTTCTTACAGAAGTCTCCTGAGCTGAGAGCTGCCAGAGTGTCTCCTTCTG
ATGTTCTGTTTGTGTTTCTTCTACTGGACAGATTGTGCTTTTCCCTAATTTAAGTCTCTCCTCAAAGG
ACAACACCTTAATGGTAAATCTCGAGAAGTCTGACCTTGGTTATGCAACGTTTAGTCCCTGGTGCT
GATTCACAGGGATGGGCTTCTGGTTGAACGTTGGCATGTTCACTGGGAGAAATTGAGAAAATGTCTCTTA
TTTATATGTTCAACAAGGAAGAAGAACTCTCAGTTCTGCAATATTGTGGGTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



ACCN:	NM_134213
Insert Size:	1104 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_134213.1, NP_598974.1</u>
RefSeq Size:	1104 bp
RefSeq ORF:	1104 bp
Locus ID:	171247
Cytogenetics:	13 A3.1