

Product datasheet for MC212400

Vmn1r205 (NM_134217) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Vmn1r205 (NM_134217) Mouse Untagged Clone
Tag: Tag Free
Symbol: Vmn1r205
Synonyms: V1rh8
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC212400 representing NM_134217
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCAAGTTTTCTACCTAGGAAGGAAGAAATGGTTTTGAATTATTAATAAAAATAATTTTTCCCTTTA
TGACTATGTTGGCACTCTGGGAACATGTCTGTTCTGTGAATTATATGTTTCAGTTGCTGGGAAGCCC
TGAGAAGAAACCCATATACCTTATTCTTATCCACTTGGCTTTTACAAACATCCTAATACTTCTTGCAAAA
GGATTGCCAAGACTATTTAGTTTTGGTTTGAGAACTTCTAGATGATATAGGGTGAAGATCATTG
TTTATCTGGAGAGGGTGGCCCGTGGCCTCTCCATCTGCACCAGCAGTCTCCTCACTGTGGTCCAGGCCAT
CATTATCAGTCCCAGAGCATCTGGTGGAGGAGGCTCAAGCCAAAGTCTGCATGGCACATTCTTCCATTC
TTCTCATTCTTTGGATACTCAATGGTTAATCAGTACGAACCTAATCCATTTTCATTAGAAGTACAAGTC
TGAATATATCACAAATATAAATGATGAGAATTATTGCTATTTTACACTACAAAATCAGGATATAAAATG
GATTTTTTAATCCTCATGGTCTGAGAGATACAGTGTTCAGGGTGCTATGGGAGGCCAGTGGCTAC
ATGGTATTTCTTCCACAAGCATCACCAGCATGTCTCTATCTTCAACTCCAATCTCCTCTACAGAA
CTCCCCGTGAGCTGAGAGCTGCTCAAAGTGTCTCCTTCTGATGCTGTGTTTTGTTTTCTTACTGGAC
AGATTGTGCCTTTTCTATAAATTATAAACATCTCTTCAAGGGACTACCCTGATGATAAAATACTCGAGAA
TTTCTGACCCCTGGTTATGCAACTTTTAGTCCCTTTGTTCTGATTACAGGGATGGGCTTCTGAATAAAT
TTTGGCATGCTCAGTGGGAGAAATTGAGAAAATATCTCTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_134217
Insert Size: 951 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_134217.1</u> , <u>NP_598978.1</u>
RefSeq Size:	951 bp
RefSeq ORF:	951 bp
Locus ID:	171251
Cytogenetics:	13 A3.1