

Product datasheet for MC212408

Vmn1r193 (NM_134225) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Vmn1r193 (NM_134225) Mouse Untagged Clone
Tag: Tag Free
Symbol: Vmn1r193
Synonyms: V1ri9
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC212408 representing NM_134225
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGTCATGTTTACTTCCTTGTGTTTCTGAAAAGGTTTTCTAGTCTGAGCATGAAGATGATTTGGAGTG
 ACCTCATCGAGGGGACAATTTCTTTCACTTGTGGACTTGGAGTTTTAGGAAACAAAATCTTATTTGT
 AAGACATTTGTATGCAGTCATCATGGGTCCTGAGAACAAAATTATAGATGTCATTCTCATCCACTGGCT
 TTTGTAACACAATCATTATTTATTGCATAGGGGTCAGAAACATAGCCACAGTTTTTTACATCAGAACT
 TCCTGGGTGATGTTGGTTGTA AAACTATAATTTATCTAGAAAAGGTTGCCCGTGGCCTCTCCATCTGCAC
 CACCTGTCTCCTCAGCGTGGTCCAGGCTGTCCACATCAGTCCCAGGACCACACTGTGGAGAAAGCTCAAA
 CCACAGACCACATGGCATGTTCTTGTCTTTCTCATCCTCTTTGGATATTTAATCCCTGATAAGCTCCA
 ACTTGCTTCACTATATCACAGCAGGCAGTAGCATGAACAGGCTGTAGTTGGGATGTTTACTGGGTATTG
 CTATATGCTGCCATCCAGGCACACAGTTAAGTGGCTGTTCTTTCTTTTTCATGGCTCTTCGTGACGTCATT
 TTCCAGGGTCTCATGGGCTGCAGCAGTGGTCTATGTCTCTCCGTCTGTACAAACATCATGTCGGTGTCC
 TCTATCTTTACAGCTCAGGGTCTGAAACAATTCAGGCCAGAAATCAGAGCTACACAGAGGGTTCTCAC
 TCTCATGACCTGTTTCCTTTTCTTTTATTGGCAGATTTCAATTTCTCCTTATACACAGGTTCCACTGTG
 ACACGTGACTCCACAATACTAAATATTAAGCATTTTTAGTACTTAGTTATGCTGGTCTCAGCCCTTTTA
 TCCTGATCATCAAGGATATTTGTGTTCTAAGACCTGCTGTGTCCTTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_134225
Insert Size: 960 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_134225.1</u> , <u>NP_598986.1</u>
RefSeq Size:	960 bp
RefSeq ORF:	960 bp
Locus ID:	171259
Cytogenetics:	13 A3.1