

Product datasheet for MC212409

Vmn1r89 (NM_134226) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Vmn1r89 (NM_134226) Mouse Untagged Clone
Tag: Tag Free
Symbol: Vmn1r89
Synonyms: V1RE8; V1rj2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC212409 representing NM_134226
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGTTTTCAAGTGACACATTCTCCAGATCTTTTCATATCCAGTTTTGTCTTGGTGTGATAGGGAAC
 CATCACTGTTAATGTTGTATGTATACACTTTCTCTTCAAGCCTCATTTTAAAAAGTTGGTAGATTCAGT
 TTTTCATGCACCTGACAATAGTTAATATGTTGACAATCGCATTACATTGATAAAAGATATCATGTTATCC
 TTTGGAGTACCAATTTCTGGATGATATCGGTTGTAAGACAGTTTTGTTTTTCATTGAGAGTCAGCCGGG
 GTCTGTCCATCTGCACCCTCTGTTCTAAGCAGTTCGAAGTCATCACCATCACTCCCAGTAACTCTAA
 GTGGGCCCTCCTGGCTCAAGCCTAGACTCTCTACATGGATTTTTCTTCTTACTTTGCTCCTGGCTTATT
 AACCTAGTCATCTACGGGTATATGCTTGACATGGTAATAGCCAAAACCAATAAAACCCATGTTGGCAATG
 GATATTCAGATGGTTACTGTCAAAACAAACACTTTGGGAACCACAATTCAGGGTCATTTTTAGCCTTCT
 CCTCACTACGATCTCTCTATGTGGCTGTCATGGTGTGGCCAGTCTCTACATGGTGATTTTTCTCTAC
 AGACTACAAGAGATCCCAGCATCTCCACAGTAGAAGCCTTCCCTCCAGTCATCTCCTGAGCGCAAAG
 CCACTCATAACATCCTGTCTCTGGTGAGCTGCTTTGTGTTAATTTATTGCTTGAACAACCATCACCCCT
 TTATGGCTTTTATACAAAAGAAAAAATCCAAGATTGGAGTCAATTAATGCAATTTTAAACAACATGCTAC
 CCAACTGTCTGTCTTTTTCTCATGAAGAATAAAAGTTATTTTGAATTCATTCTTCATTTTCGG
 TACTGAGGATGACATGTCATCAAAGTTCAGTGCATGGCTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_134226
Insert Size: 951 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_134226.1 , NP_598987.1
RefSeq Size:	951 bp
RefSeq ORF:	951 bp
Locus ID:	171260
Cytogenetics:	7 A1