

Product datasheet for **MC212410**

Vmn1r87 (NM_134227) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Vmn1r87 (NM_134227) Mouse Untagged Clone
Tag: Tag Free
Symbol: Vmn1r87
Synonyms: V1rk1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC212410 representing NM_134227
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGGAACTTGCTCCTTTCATACTATACATGTACATCTTTTTAGTTCAGCCTCATCAGAAGAAGCCCG
TAGATATGATCTTCACACATCTAATGCTTGTCAATGTTTTGAGCATCATGTTCAAGCTGATGCCAGATGT
CATGAGATCCTTTCAGTCAGGATCCTTTCCATGATGTTGTATGTAAGGCAGTTTTGTATGCATACAGT
GTTACTAGGGGCTTTCCCTCTGTACTATCTTGCTACTAAGTGTATTTCAAGCCATCACTCTCAGTTCTA
AACATTCCAAGTGGGCATGGCTTAAATCCAAGCTTGAGTCATGCATTTTTCCCTTATTCCTCATCATGTG
GATCATCAACACATTTCTCTATCTTCCCATGGTTGAAAACATAAAGGGAAAAATCAACTTCACTTTTGTG
GATTCAAAGTATTCAGATATACTGCCAAGTAGCCAGGTTCCGCCACCACACCAAGTCACTTGTAA
CTGCACTAATAATTAGAGACGTCCTGTTTGTGTTTCTCATGATATGGTCCAGCCTCTACATGGTGAGCCT
CCTCTTCAGACACAATAGGAGTACACGGCATGTCCACAGTTTCAGTGTCTCTCCATGCATCTTCTGAA
AAGAGAGCCACCCACAGCATCTTCTGCTTGGCTTTTTTCATGTTTTTTTATTGTTCAAATACCTTTG
TCACCTTCTATTCACCTTACCAACCAAGAGCAACCCAGTATTGGATCTGATTGGTGGAGTTTTATCATC
AGGCTACCCAAGTGTCTGCCCTTATGTTCTGTTGAACAATAGGAACTTATTTCCAAATTTATTTCTTCA
TTTTCAAATTTAAATGTACCTTTTTTGAAGAGGATGCCACTGCT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_134227
Insert Size: 888 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_134227.1</u> , <u>NP_598988.1</u>
RefSeq Size:	888 bp
RefSeq ORF:	888 bp
Locus ID:	171261
Cytogenetics:	7 A1